SOLUTIONS FOR SAMPLE PREPARATION







Founded in 1995, SiliCycle is specialized in the development, manufacturing and commercialization of high value silica gels and specialty products for chromatography, purification and synthesis.



Solutions for Sample Preparation

Silia <i>Prep</i> ™ & Silia <i>PrepX</i> ™ SPE Cartridges and Well Plates	2
Product Selection Guide by Technical Characteristics	3
Typical Applications & Experimental Procedures	4
Product Selection Guide by Manufacturer	7
Ordering Information	8
Applications	10
Suggested Protocols for Various Analytes & Matrices	32
Silia <i>Prep</i> ™ SPE Accessories	40
Silia <i>Prep</i> ™ Micro-SPE Tips	46
Silia <i>Fast</i> ™ for Sample Preparation & Pesticide Analysis	50
Contact Us	52





SPE Cartridges and Well Plates

- Wide variety of sorbents
- Tight particle size distribution
- Very good packing (no fines)
- High recovery and yield

Silica-based Silia*Prep*™ and Polymeric Silia*PrepX*™

Solid-phase extraction (SPE) is designed for rapid sample preparation and purification prior to chromatographic analysis.

Our SiliaPrep (silica-based) and SiliaPrepX (polymeric) families of SPE cartridges and well plates have been created to cover the entire spectrum of solid-phase extraction. This complete range of sorbents allows treatment of most common matrices:

- · human and animal biological fluids
- · petrochemical residues
- · food and beverage

· waste waters


· toxicological residues

Silia Prep and Silia PrepX products are made using state-of-the-art technology, giving you the highest quality and the best lot-to-lot reproducibility. These advanced sorbents are providing you a clean extract, reducing ion suppression and increasing selectivity for LC/MS/MS applications.

All our ultrapure SiliaFlash silica gels and functionalized SiliaBond silica gels are available in SPE formats. Just tell us what you need!

Cartridge sizes

We can provide a complete range of SPE cartridge lengths and diameters.



Bigger sizes (70 mL, 150 mL & 276 mL) are also available under SiliaSep OT branding.

Tips for your method development

	Tips for Your Method Development		
Sorbent Type	Silica-Based (SiliaPrep)	Polymeric (SiliaPrepX)	
Sorbent Capacity	Load up to 5 % of bed weight: 100 mg of silica-based sorbent will retain up to 5 mg of sample	Load up to 10 % of bed weight: 100 mg of polymeric sorbent will retain up to 10 mg of sample	
Not enough sorbent: ANALYTE LOSS ► low recovery and reproducibility Too much sorbent: MORE EXPENSIVE ► more solvent used, taller SPE cartridges Concentrated samples: double the bed weight to avoid analyte loss			



Product Selection Guide by Technical Characteristics

	Product Sel	ection Guide b	y Technica	al Characteristic	s (typical va	alues)		
SiliaPrep SiliaPrepX	Sorbent Number	Particle Size	Pore Size	Surface Area	Carbon Load	Endcapping	lonic Capacity	pH Stability
Silica-Based Non Polar P	hases							
Silia <i>Prep</i> C18 Plus	SPE-R00830B-xxx	40 - 63 μm	60 Å	500 m²/g	17 %	Proprietary	-	2 - 10
SiliaPrep C18 nec	SPE-R35530B-xxx	40 - 63 μm	60 Å	500 m ² /g	17 %	No	-	2 - 10
Silia <i>Prep</i> C18 WPD	SPE-R33229G-xxx	37 - 55 μm	125 Å	300 m²/g	13 %	Yes	-	2 - 10
Silia <i>Prep</i> C8	SPE-R31030B-xxx	40 - 63 μm	60 Å	500 m²/g	12 %	Yes	-	2 - 10
SiliaPrep C8 nec	SPE-R31130B-xxx	40 - 63 μm	60 Å	500 m²/g	12 %	No	-	2 - 10
Silia <i>Prep</i> Phenyl (<i>PH</i>)	SPE-R34030B-xxx	40 - 63 μm	60 Å	500 m²/g	9 %	Yes	-	2 - 10
Silia <i>Prep</i> PFP	SPE-R67530B-xxx	40 - 63 μm	60 Å	500 m²/g	11 %	Yes	-	2 - 10
Silica-Based Polar Phase	es							
Silia <i>Prep</i> Cyano (<i>CN</i>)	SPE-R38030B-xxx	40 - 63 μm	60 Å	500 m²/g	7 %	Yes	-	2 - 10
Silia <i>Prep</i> Diol <i>nec</i>	SPE-R35030B-xxx	40 - 63 μm	60 Å	500 m²/g	8 %	No	-	2 - 10
Silia <i>Prep</i> Florisil	SPE-AUT-0014-xxx	40 - 75 μm	100 Å	250 m²/g	-	-	-	3 - 8
Silia <i>Prep</i> Florisil LP	SPE-AUT-0014LP-xxx	75 - 150 μm	100 Å	250 m²/g	-	-	-	3 - 8
Silia <i>Prep</i> Florisil PR	SPE-AUT-0015-xxx	150 - 250 μm	-	200 m²/g	-	-	-	3 - 8
Silia <i>Prep</i> Silica	SPE-R10030B-xxx	40 - 63 μm	60 Å	500 m ² /g	-	-	-	2 - 9
Silia <i>Prep</i> Silica WPD	SPE-R10029G-xxx	37 - 55 μm	125 Å	300 m²/g	-	-	-	2 - 9
Silia <i>Prep</i> Acidic Alumina	SPE-AUT-0053-xxx	75 - 150 μm	70 Å	150 - 320 m²/g	-	-	-	3 - 8
Silia <i>Prep</i> Neutral Alumina	SPE-AUT-0054-xxx	75 - 150 μm	70 Å	150 - 320 m²/g	-	-	-	3 - 8
Silia <i>Prep</i> Basic Alumina	SPE-AUT-0055-xxx	75 - 150 μm	70 Å	150 - 320 m²/g	-	-	-	3 - 8
Silica-Based Ion Exchanç	ge Phases							
SiliaPrep SAX nec	SPE-R66530B-xxx	40 - 63 μm	60 Å	500 m²/g	10 %	No	0.90 meq/g	2 - 10
SiliaPrep SAX-2 nec	SPE-R66430B-xxx	40 - 63 μm	60 Å	500 m²/g	9 %	No	0.71 mmol/g	2 - 10
Silia <i>Prep</i> Carbonate	SPE-R66030B-xxx	40 - 63 μm	60 Å	500 m²/g	6 %	Yes	0.46 mmol/g	2 - 10
Silia <i>Prep</i> Amine (WAX)	SPE-R52030B-xxx	40 - 63 μm	60 Å	500 m²/g	7 %	Yes	1.2 mmol/g	2 - 10
Silia <i>Prep</i> SCX	SPE-R60530B-xxx	40 - 63 μm	60 Å	500 m²/g	9 %	Yes	0.54 meq/g	2 - 10
Silia <i>Prep</i> SCX-2	SPE-R51230B-xxx	40 - 63 μm	60 Å	500 m²/g	5 %	Yes	0.63 meq/g	2 - 10
Silia <i>Prep</i> WCX	SPE-R70030B-xxx	40 - 63 μm	60 Å	500 m²/g	7 %	Yes	0.92 mmol/g	2 - 10
Specialty Phases								•
Silia <i>Prep</i> PCB	SP2-R00650030B-xxx	40 - 63 μm	60 Å	500 m²/g	3 %	Proprietary	-	2 - 10
Silia <i>Prep</i> CleanDRUG	SPEC-R651230B-xxx	40 - 63 μm	60 Å	500 m²/g	9 %	Proprietary	-	2 - 10
Silia <i>Prep</i> CleanENVI	SPEC-R31930B-xxx	40 - 63 μm	60 Å	500 m²/g	19 %	Proprietary	-	2 - 10
Silia <i>Prep</i> PAH	SP2-R0610030B-xxx	40 - 63 μm	60 Å	500 m²/g	13 %	Proprietary	-	2 - 10
Polymeric Phases								
Silia <i>PrepX</i> DVB	SPE-P0001-xxx	85 μm	60 Å	1,000 m²/g	90 %	-	-	1 - 14
Silia <i>PrepX</i> HLB	SPE-P0002-xxx	40 μm	110 Å	850 m²/g	88 %	-	-	1 - 14
Silia <i>PrepX</i> SCX	SPE-P0005-xxx	85 μm	60 Å	800 m²/g	80 %	-	0.80 meq/g	1 - 14
Silia <i>PrepX</i> SAX	SPE-P0010-xxx	85 μm	60 Å	900 m²/g	85 %	-	0.20 meq/g	1 - 14
Silia <i>PrepX</i> WCX	SPE-P0015-xxx	85 μm	60 Å	800 m²/g	85 %	-	0.70 meq/g	1 - 14
Silia <i>PrepX</i> WAX	SPE-P0020-xxx	85 μm	60 Å	800 m²/g	86 %	-	0.50 meg/g	1 - 14



Reversed and Normal Phases - Typical Applications

The table below will help you select the right media to purify your compounds of interest, in either reversed-phase or normal phase.

SPE Cartridges & Well Plates Portfolio (Reversed and Normal Phases)				
Mode SiliaPrep Phases		Applications		
	SiliaPrep C18 (Plus, WPD Widepore, nec)	For organic compounds from water, drugs & metabolites from fluids		
Reversed-Phases: non-polar sorbents	Silia <i>Prep</i> C8 (endcapped & nec)	For extremely non-polar and large compounds (vitamin D, oils)		
non polar consents	SiliaPrep Phenyl (PH) & Pentafluorophenyl (PFP)	For aromatic compounds, complex natural products		
Polymeric Reversed-Phases	SiliaPrepX HLB & DVB	For drugs & metabolites from biological fluids, API from tablets and cream		
	Silia <i>Prep</i> Cyano (<i>CN</i>)	For acidic, basic and neutral compounds from aqueous solutions		
	SiliaPrep Diol nec	For polar compounds from non-polar solvents, structural isomers		
Normal Phases: polar sorbents	SiliaPrep Florisil & Florisil PR (Pesticide Residues)	For chlorinated pesticides, PCB's and polysaccharides		
	SiliaPrep Silica & Silica WPD (Widepore)	For various compounds from non-polar solvents, structural isomers		
	Silia <i>Prep</i> Alumina (<i>Acidic, Neutral & Basic</i>)	For aromatic compounds and aliphatic amines		

Experimental Procedures

Generic protocols are presented below, for reversed-phase and normal phase SPE, to help you develop your method depending on the sorbent used, the sample matrix and the analyte properties.

These are only convenient starting points for method development. Further optimization may be required to tailor the method to the application needs.

Reversed-Phases

Extraction of neutral, acidic & basic organic compounds

Extraction of neutral, acidic & basic organic compounds			
CONDITIONNING STEP	1 x CV ⁽¹⁾ of Methanol		
EQUILIBRATION STEP	1 x CV of water		
LOADING STEP	Aqueous sample, pH adjusted 2 units above pK _a (<i>bases</i>) or below pKa (<i>acids</i>)		
WASHING STEP	1 x CV of 5 % Methanol ⁽²⁾ in water		
ELUTION STEP	1 x CV of Methanol		

Normal Phases

Extraction of compounds from non-polar solvents

	Extraction of compounds from non-polar solvents
CONDITIONNING STEP	1 x CV of Isopropanol
EQUILIBRATION STEP	1 x CV of Hexane (or other low polar solvent)
LOADING STEP	Sample diluted in Hexane
WASHING STEP	1 x CV of 5 % Isopropanol in Hexane
ELUTION STEP	1 x CV of 50 - 95 % Isopropanol in Hexane

Notes:



⁽¹⁾ Abbreviation used: CV = Column Volume

⁽²⁾ For polymeric sorbents used in reversed-phase, you can add up to 60 % Methanol in water during the washing step, if your application requires it.

Ion Exchange Phases - Typical Applications

The table below will help you select the right media according to the pK_a of your analyte.

	SPE Cartridges & Well Plates Portfolio (Ion Exchange Phases)				
Mode	SiliaPrep Phases Applications				
	SiliaPrep SAX & SAX-2 (TMA Chloride & Acetate) nec	For weakly acidic molecules ($pK_a 3 - 5$)			
Ion Exchange Phases: ionic sorbents	SiliaPrep Carbonate	For scavenging of TFA, extraction of acids (boronic acids & acidic phenols)			
	SiliaPrep Amine (WAX)	For strongly acidic molecules ($pK_a < 3$), structural isomers, saccharides			
	SiliaPrep SCX & SCX-2 (Tosic & Propylsulfonic Acids)	For weakly basic molecules (pK _a 7 - 9), catch & release of amines			
	SiliaPrep WCX (Carboxylic Acid)	For strongly basic compounds ($pK_a > 9$)			
Polymeric Ion SiliaPrepX SAX & WAX		For acidic compounds & metabolites, highly stable in organic solvents			
Exchange Phases	SiliaPrepX SCX & WCX	For basic compounds, highly stable in organic solvents			

Experimental Procedures

Strong Anion Exchangers (SAX)

Extraction of weak acids (pK _a 3 - 5 ⁽³⁾)		
CONDITIONNING STEP	1 x CV of Methanol	
LOADING STEP	Aqueous sample, pH adjusted at 7.0 - 8.0	
WASHING STEP	1 x CV of Methanol (elution of basic & neutral compounds)	
ELUTION STEP	1 x CV of 2 - 5 % HCO₂H in Methanol (elution of weak acidic compounds)	

Strong Cation Exchangers (SCX)

Extraction of weak bases (pK _a 7 - 9)			
CONDITIONNING STEP	1 x CV of Methanol		
EQUILIBRATION STEP	1 x CV of water		
LOADING STEP	Aqueous sample, pH adjusted at 3.0 - 4.0		
WASHING STEP 1	1 x CV of water		
WASHING STEP 2	1 x CV of Methanol (elution of acidic & neutral compounds)		
ELUTION STEP	1 x CV of 2 - 5 % NH₄OH ⁽⁴⁾ in Methanol (elution of weak basic compounds)		

Weak Anion Exchangers (WAX)

Extraction of strong acids (pK _a < 3)		
CONDITIONNING STEP	1 x CV of Methanol	
EQUILIBRATION STEP	1 x CV of water	
LOADING STEP	Aqueous sample, pH adjusted at 4.0 - 5.0	
WASHING STEP 1	1 x CV of water	
WASHING STEP 2	1 x CV of Methanol (elution of basic & neutral compounds)	
ELUTION STEP	1 x CV of 2 - 5 % NH₄OH ⁽⁴⁾ in Methanol (elution of strong acidic compounds)	

Weak Cation Exchangers (WCX)

Extraction of strong bases ($pK_a > 9$)		
CONDITIONNING STEP	1 x CV of Methanol	
EQUILIBRATION STEP	1 x CV of water	
LOADING STEP	Aqueous sample, pH adjusted at 8.0	
WASHING STEP 1	1 x CV of water	
WASHING STEP 2	1 x CV of Methanol (elution of acidic & neutral compounds)	
ELUTION STEP	1 x CV of 2 - 5 % HCO ₂ H in Methanol (elution of strong basic compounds)	

Notes:



⁽³⁾ For extraction of Phenol (pK, 10), we recommend using a polymeric phase (SiliaPrepX SAX) and load the sample with a pH adjusted to 12.

⁽⁴⁾ For silica-based sorbents, NH₄OH can be too aggressive. You can use NH₂ (7M) in Methanol to avoid degrading the phase.

Specialty Phases & Metal Scavengers - Typical Applications

The table below presents our specialty phases, to remove specific compounds from your samples.

4	SPE Cartridges & Well Plates Portfolio (Specialty Phases & Metal Scavengers)			
Mode	SiliaPrep Phases	Applications		
	SiliaPrep PCB	For extraction of PCB's from waste oil (hexane extract)		
Cussialty Dhases	Silia <i>Prep</i> CleanDRUG	For drugs of abuse applications		
Specialty Phases	SiliaPrep CleanENVI	For PAH's, PCB's, herbicides and pesticides from waste waters		
	SiliaPrep PAH	For PAH's from waste waters		
Matal Scavandare		For lowering the residual metal concentration of various metal complexes (<i>Pd</i> , <i>Pt</i> , <i>Rh</i> , <i>Ru</i> , <i>Ni</i> , <i>Sn</i> , <i>etc</i>) to single digit ppm		

Experimental Procedures

The procedures below are only convenient starting points for method development. Further optimization may be required to tailor the method to your application needs.

Specialty Phases

Extraction of PCBs, drugs and PAHs

PCBs from waste oil with SiliaPrep PCB			
CONDITIONNING STEP	1 x CV ⁽¹⁾ of Hexane		
LOADING STEP	Diluted sample (with Hexane)		
ELUTION STEP	1 x CV of Hexane		

Metal Scavengers

Catch of the metal & release of your analyte

Catch and release of the analyte				
EQUILIBRATION STEP	1 x CV of sample solvent			
LOADING STEP(2)	Diluted sample			
RINSING STEP	1 x CV of sample solvent			

Drugs of abuse with SiliaPrep CleanDRUG			
CONDITIONNING STEP	1 x CV of Methanol		
EQUILIBRATION STEP	1 x CV of water (buffered at pH 6.0)		
LOADING STEP	Aqueous sample (buffered at pH 6.0)		
WASHING STEP	1 x CV of water then 1 x CV of Methanol		
ELUTION STEP	1 x CV of Isopropanol:NH ₄ OH (90:10)		

Environmental samples with SiliaPrep CleanENVI & PAH			
CONDITIONNING STEP	$1 \times \text{CV}$ of Dichloromethane then $1 \times \text{CV}$ of Methanol		
EQUILIBRATION STEP	1 x CV of water		
LOADING STEP	Aqueous sample		
WASHING STEP	1 x CV of 5 - 50 % Methanol in water		
ELUTION STEP	1 x CV of Dichloromethane		

Notes



⁽¹⁾ Abbreviation used: CV = Column Volume

⁽²⁾ Non retentive SPE (Catch & Release): analyte won't retain on the sorbent and will elute directly during loading and rinsing steps. Scavenged compounds will remain in the SPE cartridge.

Product Selection Guide by Manufacturer

The table below will help you find equivalences to products of well-known SPE manufacturers.

Product Selection Guide by Manufacturer					
SiliCycle	Waters	Phenomenex	Agilent	Biotage	Macherey-Nagel
Silia <i>Prep</i> C18 Plus	Sep-Pak® tC18	Strata® C18-E	Bond Elut® C18	Isolute® C18 (EC)	Chromabond® C18 ec
Silia <i>Prep</i> C18 nec		Strata® C18-U	Bond Elut® C18 OH	Isolute® C18	Chromabond® C18
Silia <i>Prep</i> C18 WPD	Sep-Pak® C18	Strata® C18-T	Bond Elut® C18 EWP	Isolute® MFC18	Chromabond® C18 ec f
Silia <i>Prep</i> C8	Sep-Pak® C8	Strata® C8	Bond Elut® C8	Isolute® C8 (EC)	
Silia <i>Prep</i> C8 nec				Isolute® C8	Chromabond® C8
Silia <i>Prep</i> Phenyl (<i>PH</i>)		Strata® Phenyl	Bond Elut® PH	Isolute® PH	Chromabond® C ₆ H ₅
Silia <i>Prep</i> PFP					
Silia <i>Prep</i> Cyano (<i>CN</i>)	Sep-Pak® Cyanopropyl	Strata® CN	Bond Elut® Cyano (CN)	Isolute® CN	Chromabond® CN
Silia <i>Prep</i> Diol <i>nec</i>	Sep-Pak® Diol		Bond Elut® Diol (20H)	Isolute® DIOL	Chromabond® OH (<i>Diol</i>)
Silia <i>Prep</i> Silica		Strata® Silica (Si-1)	Bond Elut® SI	Isolute® SI	Chromabond® SiOH
Silia <i>Prep</i> Silica WPD	Sep-Pak® Silica				
Silia <i>Prep</i> Florisil LP & Florisil PR	Sep-Pak® Florisil®	Strata® FL-PR (Florisil®)	Bond Elut® Florisil	Isolute® FL	Chromabond® Florisil®
Silia <i>Prep</i> Alumina (<i>Acidic, Neutral, Basic</i>)	Sep-Pak [®] Alumina (A, N, B)	Strata® Alumina-N (<i>AL-N</i>)	Bond Elut® Alumina (-A, -N, -B)	Isolute® ALUMINA (AL-A, AL-N & AL-B)	Chromabond® Alox (A, N, B)
SiliaPrep SAX nec (TMA Chloride)	Sep-Pak® Accell™ Plus QMA	Strata® SAX	Bond Elut® SAX	Isolute® SAX	Chromabond® SB
SiliaPrep SAX-2 nec (TMA Acetate)				Isolute® PE-AX	
Silia <i>Prep</i> Carbonate	Accell Plus QMA Carbonate Plus Light			Isolute® Si-Carbonate (Si-TMA-CO ₃)	
Silia <i>Prep</i> Amine (<i>WAX</i>)	Sep-Pak® Amino	Strata® NH ₂	Bond Elut® NH2	Isolute® NH ₂	Chromabond® NH ₂
Silia <i>Prep</i> Tosic Acid (SCX)		Strata® SCX	Bond Elut® SCX	Isolute® SCX-3	Chromabond® SA
SiliaPrep SCX-2 (Propylsulfonic Acid)			Bond Elut® PRS	Isolute® SCX-2	Chromabond® PSA
SiliaPrep WCX (Carboxylic Acid)	Sep-Pak® Accell™ Plus CM	Strata® WCX	Bond Elut® CBA	Isolute® CBA	Chromabond® PCA
Silia <i>Prep</i> PCB			Bond Elut® PCB		Chromabond® SA/SiOH
Silia <i>Prep</i> CleanDRUG		Strata® Screen-C	Bond Elut® Certify	Isolute® HCX	Chromabond® Drug
Silia <i>Prep</i> CleanENVI & Silia <i>Prep</i> PAH		Strata® PAH	EnvirElut®	Isolute® PAH	Chromabond® C18 PAH
Silia <i>PrepX</i> HLB	Oasis® HLB	Strata®-X	Bond Elut® NEXUS		Chromabond® HLB
Silia <i>PrepX</i> DVB			Bond Elut® ENV	Isolute® 101	Chromabond® HR-X
Silia <i>PrepX</i> SAX	Oasis® MAX	Strata®-X-A	Bond Elut® Plexa PAX	Evolute® Express AX	Chromabond® HR-XA
Silia <i>PrepX</i> WAX	Oasis® WAX	Strata®-X-AW		Evolute® Express WAX	Chromabond® HR-XAW
Silia <i>PrepX</i> SCX	Oasis® MCX	Strata®-X-C	Bond Elut® Plexa PCX	Evolute® Express CX	Chromabond® HR-XC
Silia <i>PrepX</i> WCX	Oasis® WCX	Strata®-X-CW	Bond Elut® NEXUS WCX	Evolute® Express WCX	Chromabond® HR-XCW

All trademarks and registered trademarks are the property of their respective owners. SiliCycle takes no responsibility for any error or omission relating to this information.



Ordering Information

To build your own product number, just add the **Phase** to the **Format PN** Examples:

- SPE-R67530B-06P for SiliaPrep PFP, 6 mL / 500 mg cartridges
- SPE-P0002-12S for SiliaPrepX HLB, 12 mL / 500 mg cartridges
- SPEC-R31930B-06S for Silia Prep Clean ENVI, 6 mL / 1 g cartridges



Silica-based Phases

SiliaPrep Phases			
Phases	Code		
Reversed-phases			
C18 Plus	R00830B		
C18 WPD	R33229G		
C18 nec	R35530B		
C8	R31030B		
C8 nec	R31130B		
Phenyl (PH)	R34030B		
PFP	R67530B		
Normal Phases	· ·		
Cyano (CN)	R38030B		
Diol nec	R35030B		
Florisil	AUT-0014		
Florisil LP	AUT-0014LP		
Florisil PR	AUT-0015		
Silica	R10030B		
Silica WPD	R10029G		
Acidic Alumina	AUT-0053		
Neutral Alumina	AUT-0054		
Basic Alumina	AUT-0055		
Ion Exchange Ph	ases		
SAX nec	R66530B		
SAX-2 nec	R66430B		
Carbonate	R66030B		
Amine (WAX)	R52030B		
Tosic Acid (SCX)	R60530B		
SCX-2	R51230B		
WCX	R70030B		
Scavengers	į.		
Cysteine	R80530B		
DMT	R79030B		
TAAcOH	R69030B		
TAAcONa	R69230B		
Thiol	R51030B		
Thiourea	R69530B		
Imidazole	R79230B		
Triamine	R48030B		

Cilia Dran Formata				
	Silia <i>Prep</i> Formats			
Formats	Qty/Box	Format PN		
SiliaPrep SPE Cartrid	ges			
1 mL / 50 mg	100	SPE-PHASE-01B		
1 mL / 100 mg	100	SPE-PHASE-01C		
3 mL / 200 mg	50	SPE-PHASE-03G		
3 mL / 500 mg	50	SPE-PHASE-03P		
6 mL / 500 mg	50	SPE-PHASE-06P		
6 mL / 1 g	50	SPE-PHASE-06S		
6 mL / 2 g	50	SPE-PHASE-06U		
12 mL / 2 g	20	SPE-PHASE-12U		
25 mL / 5 g*	20	SPE-PHASE-20X		
70 mL / 10 g*	16	FLH-PHASE-70Y		
70 mL / 15 g*	16	FLH-PHASE-70i		
70 mL / 20 g*	16	FLH-PHASE-70Z		
150 mL / 25 g*	10	FLH-PHASE-95K		
150 mL / 50 g*	10	FLH-PHASE-95M		
150 mL / 70 g*	10	FLH-PHASE-95N		
276 mL / 100 g*	12	FLH-PHASE-276F		
SiliaPrep Large Rese	rvoir Volume SPE Cartridge	es		
10 mL / 200 mg	50	SPC-PHASE-10G		
10 mL / 500 mg	50	SPC-PHASE-10P		
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-PHASE-P		
1 g	50	SPS-PHASE-S		
SiliaPrep 96-Well Plates				
2 mL / 50 mg	1	96W-PHASE-B		
2 mL / 100 mg	2 mL / 100 mg 1			



Polymeric Phases

SiliaPrepX Phases		
Phase	Code	
DVB	P0001	
HLB	P0002	
SCX	P0005	
SAX	P0010	
WCX	P0015	
WAX	P0020	

SiliaPrepX Formats					
Formats	Qty/Box	Format PN			
Silia <i>PrepX</i> SPE Car	SiliaPrepX SPE Cartridges				
1 mL / 30 mg	100	SPE-PHASE-01AA			
3 mL / 30 mg	50	SPE-PHASE-03AA			
3 mL / 60 mg	50	SPE-PHASE-03BB			
6 mL / 100 mg	30	SPE-PHASE-06C			
6 mL / 200 mg	30	SPE-PHASE-06G			
6 mL / 500 mg	30	SPE-PHASE-06P			
12 mL / 500 mg	20	SPE-PHASE-12P			
12 mL / 1 g	20	SPE-PHASE-12S			
25 mL / 1 g*	20	SPE-PHASE-20S			
25 mL / 2 g*	20	SPE-PHASE-20U			
70 mL / 5 g*	16	FLH-PHASE-70X			
SiliaPrepX 96-Well Plates					
2 mL / 10 mg	1	96W-PHASE-1A			
2 mL / 30 mg	1	96W-PHASE-AA			
2 mL / 60 mg	1	96W-PHASE-BB			

Specialty Phases

Silia <i>Prep</i> Specialty Phases			
Phase	Code		
PCB	R00650030B		
PAH	R0610030B		
CleanDRUG	R651230B		
CleanENVI	R31930B		

SiliaPrep Specialty Formats					
Formats	Qty/Box	PCB & PAH	CleanDRUG & CleanENVI		
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SP2-PHASE-01B	SPEC-PHASE-01B		
1 mL / 100 mg	100	SP2-PHASE-01C	SPEC-PHASE-01C		
3 mL / 200 mg	50	SP2-PHASE-03G	SPEC-PHASE-03G		
3 mL / 500 mg	50	SP2-PHASE-03P	SPEC-PHASE-03P		
6 mL / 500 mg	50	SP2-PHASE-06P	SPEC-PHASE-06P		
6 mL / 1 g	50	SP2-PHASE-06S	SPEC-PHASE-06S		
6 mL / 2 g	50	SP2-PHASE-06U	SPEC-PHASE-06U		
12 mL / 2 g	20	SP2-PHASE-12U	SPEC-PHASE-12U		
25 mL / 5 g*	20	SP2-PHASE-20X	SPEC-PHASE-20X		
70 mL / 10 g*	16	FLH-PHASE-70Y			
70 mL / 15 g*	16	FLH-PHASE-70i			
70 mL / 20 g*	16	FLH-PHASE-70Z			
150 mL / 25 g*	10	FLH-PHASE-95K			
150 mL / 50 g*	10	FLH-PH/	ASE-95M		
150 mL / 70 g*	10	FLH-PHASE-95N			
276 mL / 100 g*	12	FLH-PHASE-276F			
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-PHASE-10G			
10 mL / 500 mg	50	SPC-PHASE-10P			
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-PHASE-P			
1 g	50	SPS-PHASE-S			

^{*} Commercialized under SiliaSepX OT branding.



Silia*Prep*[™] and Silia*PrepX* [™] Applications

Forensics

Opioids	11
Antidepressants	12
Pharmaceutical drugs	13
Anti-inflammatory drugs	14
Performance-enhancing drugs	15
Cannabinoids	16
Stimulants	17
Varia	18

Food

Orange & apple juices	19
Milk & honey	20
Salmon & beef	21
Potatoes	22

Environment

Pesticides in water	23
Bisphenol A in water	26
Drugs in water	27
Other analytes in water	28
Other matrices	30







Extraction of Methadone from Human Urine and Serum CARTRIDGE SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G SAMPLE PRETREATMENT 200 µL of Phosphoric Acid 2 % was added to 1 mL of urine / serum sample CONDITIONNING STEP 6 mL of Methanol **EQUILIBRATION STEP** 6 mL of Water LOADING STEP Treated sample was slowly aspirated through the cartridge **WASHING STEP** 6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge **ELUTION STEP** 2 x 3 mL of 10 % Ammonia in Methanol FURTHER TREATMENT Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS at 1 µg/mL RECOVERY Methadone in urine 90 %

Methadone in serum

95 %

Extraction of Methadone and EDDP from Human Urine			
CARTRIDGE	Silia <i>PrepX</i> HLB 1 mL / Part Number: SPE-P00	0	
SAMPLE PRETREATMENT	40 µL of internal standa 200 ng/mL in Methanol of urine sample and 20 Hydroxide 4 %) was added to 200 μL	
CONDITIONNING STEP	1 mL of Methanol		
EQUILIBRATION STEP	1 mL of Ammonium Hy	droxide 2 %	
LOADING STEP	Urine sample was slowly aspirated through the cartridge		
WASHING STEP	1 mL of Methanol / Ammonium Hydroxide 2 % (50:50) then 1 mL of Methanol / Ammonium Hydroxide 2 % (80:20)		
ELUTION STEP	1 mL of Methanol / water (80:20)		
FURTHER TREATMENT	Quantification by LDTD/MS/MS (collaboration with Phytronix)		
RECOVERY	at 1,000 ng/mL		
	Methadone	91 %	
	EDDP	85 %	

Extraction of Fen	tanyl and Norfentany	l from Urine
CARTRIDGE	Silia <i>Prep</i> CleanDRUG Part Number: SPEC-R	9
SAMPLE PRETREATMENT	200 μL of urine was added to 600 μL of Sodium Acetate in water and 40 μL of internal standard (200 ng/mL in Methanol)	
CONDITIONNING STEP	1 mL of Methanol	
EQUILIBRATION STEP	1 mL of Water and 1 mL of Sodium Acetate in Water (100 mM, pH 6.0)	
LOADING STEP	Urine sample was slowly aspirated through the cartridge	
WASHING STEP	1 mL of Water then 1 mL of Methanol	
ELUTION STEP	1 mL of Ethyl Acetate / Isopropanol / Ammonium Hydroxide (78:20:2)	
FURTHER TREATMENT	Evaporation, reconstitution and quantification by LDTD/MS/MS (collaboration with Phytronix)	
RECOVERY	at 500 ng/mL	
	Fentanyl	96 %
	Norfentanyl	98 %

Extraction of C	Codeine from Human Urine and Serum	
CARTRIDGE	Silia <i>PrepX</i> SCX 6 mL / 200 mg Part Number: SPE-P0005-06G	
SAMPLE PRETREATMENT	200 μL of Phosphoric Acid 2 % was added to 1 mL of urine / serum sample	
CONDITIONNING STEP	6 mL of Methanol	
EQUILIBRATION STEP	6 mL of water	
LOADING STEP	Treated sample was slowly aspirated through the cartridge	
WASHING STEP	6 mL of Hydrochloric Acid 0.1 N then 6 mL of Methanol, dry the cartridge	
ELUTION STEP	2 x 3 mL of 5 % Ammonia in Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS	
RECOVERY	at 1 μg/mL	
	Codeine in urine 70 %	
	Codeine in serum 92 %	





Extraction of Tricyclic Antidepressants from Serum			
CARTRIDGE	Silia <i>PrepX</i> WCX 3 mL Part Number: SPE-P00	0	
SAMPLE PRETREATMENT	250 μL of serum were diluted with 1 mL of 10 % Formic Acid in Water		
CONDITIONNING STEP	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of Water		
LOADING STEP	Treated sample was slowly aspirated through the cartridge		
WASHING STEP	1 mL of 5 % Formic Acid in Water then 1 mL Methanol, dry the cartridge		
ELUTION STEP	3 mL of 5 % Formic Acid in Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS		
RECOVERY	at 1 µg/mL		
	Doxepine	79 %	
	Imipramine	79 %	
	Amitriptyline	91 %	
	Trimipramine	98 %	

Extraction o	f Pharmaceutic	al Drugs	from Serum	
CARTRIDGE	Silia <i>PrepX</i> SCX 6 mL / 200 mg Part Number: SPE-P0005-06G			
SAMPLE PRETREATMENT	200 μL of Phosphoric Acid 2 % was added to 1 mL of serum sample) 1 mL
CONDITIONNING STEP	6 mL of Methan	ol		
EQUILIBRATION STEP	6 mL of Water			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	6 mL of Chlorhydric Acid 0.1 N, dry the cartridge			
ELUTION STEP	2 x 3 mL of Methanol (<i>acidic and neutrals</i> analytes) and 2 x 3 mL of 10 % Ammonia in Methanol (<i>basic analytes</i>)			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS			
RECOVERY	at 1 μg/mL			
	Indomethacin	33 %	Phenobarbital	108 %
	Tolmetin	73 %	Trimipramine	92 %
	Hexobarbital	80 %	Amitriptyline	94 %
	Naproxen	85 %	Imipramine	95 %
	Suprofen	108 %	Doxepin	101 %

Ropinirole & Ar	nitriptyline Detection	in Human Plasma
CARTRIDGE	Silia <i>Prep</i> CleanDRUG Part Number: SPEC-R	S .
SAMPLE PRETREATMENT	Mix 0.1 mL of plasma w Water (50:50) and 2 m	vith 0.1 mL of Methanol and L of 1 % Acetic Acid
CONDITIONNING STEP	3 mL of Methanol	
EQUILIBRATION STEP	3 mL of Water	
LOADING STEP	Plasma sample was sle cartridge	owly aspirated through the
WASHING STEP	3 mL of Water then 3 n	nL of Methanol
ELUTION STEP	3 mL of 5 % Ammonium Hydroxide in Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS	
RECOVERY	at 10 ng/mL	
	Ropinirole	94 %
	Amitriptyline	90 %

Futuration o	f Dhawmaaautic	ol Davago	from Corum	
Extraction o	f Pharmaceutic	ai Drugs	from Serum	
CARTRIDGE	Silia <i>PrepX</i> SAX 6 mL / 200 mg Part Number: SPE-P0010-06G			
SAMPLE PRETREATMENT	pH of serum sample was adjusted to basic value with Sodium Hydroxide 1 N			value
CONDITIONNING STEP	6 mL of Methan	ol		
EQUILIBRATION STEP	6 mL of Water			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	6 mL of Water, dry the cartridge			
ELUTION STEP	2 x 3 mL of Methanol (<i>basic analytes</i>) and 2 x 3 mL of Formic Acid 10 % in Methanol (<i>acidic</i> <i>analytes</i>)			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS			
RECOVERY	at 1 µg/mL			
	Nortriptyline	69 %	Imipramine	80 %
	Doxepine	72 %	Tolmetin	85 %
	Trimipramine	73 %	Naproxen	86 %
	Protriptyline	75 %	Suprofen	96 %
	Amitriptyline	78 %		





Extraction of	Tricyclic Antidepress	ants from Serum
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G	
CONDITIONNING STEP	5 mL of Methanol	
EQUILIBRATION STEP	5 mL of Water	
LOADING STEP	1 mL of serum sample (pH value adjusted with 25 μL of Phosphoric Acid) was slowly aspirated through the cartridge	
WASHING STEP	5 mL of Water, dry the cartridge	
ELUTION STEP	2 x 3 mL of Methanol	
FURTHER TREATMENT	Quantification by LC/MS	
RECOVERY	Protriptyline	80 %
	Nortriptyline	75 %
	Doxepine	91 %
	Imipramine	88 %
	Amitriptyline	88 %
	Trimipramine	88 %

Extraction of Barbiturates from Serum			
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / Part Number: SPE-P0	•	
CONDITIONNING STEP	6 mL of Methanol		
EQUILIBRATION STEP	6 mL of Water		
LOADING STEP	1 mL of serum sample through the cartridge	1 mL of serum sample was slowly aspirated through the cartridge	
WASHING STEP	6 mL of Water, dry the cartridge		
ELUTION STEP	6 x 1 mL of Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS		
RECOVERY	at 100 ng/ml		
	Phenobarbital	99 %	
	Pentobarbital	69 %	
	Hexobarbital	86 %	

Extraction of Antibacterial Drugs from Serum					
CARTRIDGE	Silia <i>PrepX</i> DVB 6 r Part Number: SPE				
SAMPLE PRETREATMENT	Mix 0.1 mL of plasr	ma with 2 mL of 1 % Acetic Acid			
CONDITIONNING STEP	6 mL of Methanol				
EQUILIBRATION STEP	6 mL of Water				
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge				
WASHING STEP	6 mL of Water, dry the cartridge				
ELUTION STEP	2 x 3 mL of Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS				
RECOVERY	Cinoxacin	100 %			
	Penicillin G	94 %			
	Penicillin V	90 %			
	Cloxacillin	88 %			

Sibutramine Detection in Human Plasma				
CARTRIDGE	Silia <i>Prep</i> CleanDRUG 3 mL / 500 mg Part Number: SPEC-R651230B-03P			
SAMPLE PRETREATMENT	Mix 0.1 mL of plasma with 2 mL of 1 % Acetic Acid			
CONDITIONNING STEP	3 mL of Methanol			
EQUILIBRATION STEP	3 mL of Water			
LOADING STEP	Plasma sample was slowly aspirated through the cartridge			
WASHING STEP	3 mL of Water then 3 mL of Methanol			
ELUTION STEP	3 mL of 5 % Ammonium Hydroxide in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS			
RECOVERY	at 5 ng/mL: 82 %			





		Extraction	of Acidic Pharmace	uticals from Se	erum
CARTRIDGE	Silia <i>PrepX</i> SAX 6 mL / 200 mg Part Number: SPE-P0010-06G				
SAMPLE PRETREATMENT	pH of serum sample was adjusted to basic value with Sodium Hydroxide 1 N				
CONDITIONNING STEP	6 mL of Methanol				
EQUILIBRATION STEP	6 mL of Water				
LOADING STEP	Treated sample was slowly aspirated through the cartridge				
WASHING STEP	6 mL of Water, then 6 mL of Sodium Hydroxide 0.1 N and 6 mL of Methanol, dry the cartridge				
ELUTION STEP	6 mL of 1 % Formic Acid in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS				
RECOVERY	at 1 µg/mL				
	Carprofen	69 %	Diclofenac	95 %	
	Ibuprofen	88 %	Fenoprofen	98 %	
	Ketoprofen	90 %	Fenoprop	104 %	
	Meclofenamic Acid	92 %	Flurbiprofen	106 %	

	E	xtraction c	of Anti-inflammatory Drug	gs From Serum
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G			
CONDITIONNING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of Water			
LOADING STEP	1 mL of serum sample (pH value adjusted with 25 μL of Phosphoric Acid) was slowly aspirated through the cartridge			
WASHING STEP	6 mL of 5 % Metanol in Water, dry the cartridge			
ELUTION STEP	2 x 3 mL of Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS			
RECOVERY	Suprofen	89 %	Naproxen	87 %
	Tolmetin	89 %	Flurbiprofen	87 %
	Sulindac	84 %	Indomethazin	85 %
	Piroxicam	86 %	Acetyl Salicylic Acid	72 %





Determination of Clenbuterol in Human Plasma				
CARTRIDGE	Silia <i>Prep</i> CleanDRUG 1 mL / 100 mg Part Number: SPEC-R651230B-01C			
SAMPLE PRETREATMENT	50 μL of internal standard (<i>Clenbuterol-d9 at 20 ng/mL in Methanol</i>) was added to 500 μL of plasma and 500 μL of Sodium Acetate (<i>100 mM</i> , <i>pH 6.0</i>)			
CONDITIONNING STEP	1 mL of Methanol			
EQUILIBRATION STEP	1 mL of Water and 1 mL of Sodium Acetate (100 mM, pH 6.0)			
LOADING STEP	Plasma sample was slowly aspirated through the cartridge			
WASHING STEP	1~mL of Water, then $1~mL$ of Acetic Acid $1~M$ and $2~x~1~mL$ of Methanol			
ELUTION STEP	1 mL of Ethyl Acetate / Isopropanol / Ammonium Hydroxide (78:20:2)			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LDTD/MS/MS (collaboration with Phytronix)			
RECOVERY	at 100 pg/mL: 94 %			

Extraction of Atenolol from Human Urine				
CARTRIDGE	Silia <i>PrepX</i> WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB			
SAMPLE PRETREATMENT	9 mL of urine was mixed with 1 mL of a solution of Atenolol in Methanol / Water (10:90)			
CONDITIONNING STEP	2 mL of Methanol			
EQUILIBRATION STEP	2 mL of Water			
LOADING STEP	1 mL of sample solution was slowly aspirated through the cartridge			
WASHING STEP	2 mL of Monopotassium Phosphate 25 mM (<i>pH 5.0</i>) then 2 mL of Methanol, dry the cartridge			
ELUTION STEP	2 mL of 2 % Formic Acid in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS			
RECOVERY	at 10 µg/mL: 90 %			

Determination of Testosterone in Human Urine			
CARTRIDGE	Mini-Silia <i>Prep</i> C18 WPD 500 mg Part Number: SPS-R33229G-P		
CONDITIONNING STEP	5 mL of Methanol		
EQUILIBRATION STEP	5 mL of Water		
LOADING STEP	2mL of urine sample was slowly aspirated through the cartridge		
WASHING STEP	5 mL of Water then 5 mL of Hexane		
ELUTION STEP	5 mL of Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using Girard-P and quantification by LC/MS/MS		
RECOVERY	at 250 ng/mL: 95 %		

Extra	Extraction of Steroids from Serum		
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G		
CONDITIONNING STEP	5 mL of Methanol		
EQUILIBRATION STEP	5 mL of Water		
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge		
WASHING STEP	5 mL of 5 % Metanol in Water, dry the cartridg	е	
ELUTION STEP	2 x 3 mL of Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution wit Acetonitrile / Water and quantification by LC/N		
RECOVERY	Methyl-6a-hydroxy-11ß-progesterone	89 %	
	Methyl-6a-hydroxy-17a-progesterone	86 %	
	Methyl-6a-hydroxy-17a-progesterone acetate	84 %	
	Hydrocortisone-21-acetate	31 %	





4		Isolation of Synthetic	Cannabinoid Metabolit	es from Urine	
CARTRIDGE	Silia <i>Prep</i> CleanDRUG 1 mL / 30 mg Part Number: SPEC-R651230B-03G				
SAMPLE PRETREATMENT	1 mL of synthetic solution (pH 6.0)	·	etabolites and deuterated i	nternal standard, then diluted	with 2 mL of a Phosphate buffer
CONDITIONNING STEP	3 mL of Methano	ol			
EQUILIBRATION STEP	3 mL of Water ar	nd 1 mL of Phosphate buffer			
LOADING STEP	Urine sample was	s slowly aspirated through the c	cartridge		
WASHING STEP	3 mL of Water then 3 mL of Phosphate buffer / Acetonitrile (80:20)				
ELUTION STEP	6 mL of Ethyl Acetate / Methanol (90:10)				
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using BSTFA and TMCS, and quantification by GC/MS				
RECOVERY	at 1,000 ng/mL	SiliaPrepX Clean DRUG	Bond Elut® Certify II	HyperSep [™] Verify AX	Clean Screen® CSTHC
	JWH-018	102 %	109 %	112 %	97 %
	JWH-122	96 %	72 %	111 %	80 %
	JWH-250 101 % 71 % 118 % 89 %				
CONCLUSION	Our Silia <i>Prep</i> Cl	eanDRUG performs as well a	s competitive products to e	extract cannabinoid metabolite	es from urine.

Source: Thesis "An Evaluation of Commercially Available Solid Phase Extraction Cartridges for the Isolation of Synthetic Cannabinoid Metabolites from urine", by Amanda Marie Forni, B.S., Ohio University, **2011**

	Detection of Δ ⁹ -Tetrahydrocannabinol in Human Plasma		
CARTRIDGE	SiliaPrep CleanENVI 3 mL / 500 mg Part Number: SPEC-R31930B-03P		
SAMPLE PRETREATMENT	250 μL of plasma was added to 1 mL Phosphate buffer (0.1 M, pH 6.0)		
CONDITIONNING STEP	3 mL of Methanol, then 3 mL of Hydrochloric Acid 1 M and 3 mL of Water		
EQUILIBRATION STEP	5 mL of Water		
LOADING STEP	Plasma sample was slowly aspirated through the cartridge		
WASHING STEP	2 mL of Water, then 1 mL of Acetic Acid 1 M and 2 mL of 20 % Methanol in Water		
ELUTION STEP	5 mL of Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using Dansyl Chloride, liquid-liquid extraction, centrifugation, evaporation under Nitrogen, reconstitution with Formic Acid / Acetone and quantification by LC/MS		
RECOVERY	at 2 ng/mL		
	THC 80 %		
	<u>THC-COOH</u> 99 %		
	THC-OH 92 %		





Extraction of Phencyclidine (PCP) from Human Urine				
CARTRIDGE	Silia <i>PrepX</i> HLB 1 mL / 30 mg Part Number: SPE-P0002-01AA			
SAMPLE PRETREATMENT	40 μL of internal standard (<i>PCP-d5 at 200 ng/mL in Methanol</i>) was added to 200 μL of urine sample and 200 μL of Ammonium Hydroxide 4 %			
CONDITIONNING STEP	1 mL of Methanol			
EQUILIBRATION STEP	1 mL of Ammonium Hydroxide 2 %			
LOADING STEP	Urine sample was slowly aspirated through the cartridge			
WASHING STEP	1 mL of Methanol / Ammonium Hydroxide 2 % (50:50) then 1 mL of Methanol / Ammonium Hydroxide 2 % (80:20)			
ELUTION STEP	1 mL of Methanol / Hydrochloric Acid 0.02 N (80:20)			
FURTHER TREATMENT	Quantification by LDTD/MS/MS (collaboration with Phytronix)			
RECOVERY	at 25 ng/mL: 99 %			

Drugs of Abuse Determination in Human Urine					
CARTRIDGE	Silia <i>Prep</i> CleanDRUG 3 mL / 200 mg Part Number: SPEC-R651230B-03G				
SAMPLE PRETREATMENT	0.5 mL of urine sample was mixed with 2.5 mL Sulfuric Acid 0.1 M				
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of Sulfuric Acid (0.1 M			
LOADING STEP	2 mL of urine sample was slowly aspirated through the cartridge				
WASHING STEP	3 mL of Phosphate buffer (pH 7.0), then 3 mL of Sulfuric Acid 0.1 M and 3 mL of Methanol				
ELUTION STEP	2 x 3 mL of Ammonium Hydroxide (5 % in Methanol)				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Water / Methanol and quantification by LC/MS				
RECOVERY	at 25 ng/mL				
	MDMA	92 %			
	MDEA	89 %			
	Amphetamine	82 %			

	An	nphetamine Quantif	ication in Human Urin	e	
CARTRIDGE	Silia <i>PrepX</i> HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
SAMPLE PRETREATMENT	100 μL of TFA was added to 10	mL of urine			
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of Water	3 mL of Water			
LOADING STEP	1 mL of urine sample was slowly	aspirated through the ca	artridge		
WASHING STEP	3 mL of (5:95) Methanol / Water with 2 % Ammonium Hydroxide; then 3 mL of (20:80) Methanol / Water with 2 % Ammonium Hydroxide and 1 mL of (80:20) Methanol / Water				
ELUTION STEP	3 mL of Methanol then 3 mL of 2 % Formic Acid in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, red	constitution with Metha	nol / Water (70:30) and o	uantification by LC/MS	3
RECOVERY	at 100 ng/mL	Silia <i>PrepX</i> HLB	Bond Elut® Plexa	Oasis® HLB	Strata [™] -X
	Amphetamine	91 %	88 %	75 %	87 %
	MDA	86 %	86 %	91 %	98 %
	MDEA	95 %	97 %	90 %	101 %
	MDMA	92 %	94 %	91 %	101 %
	Methamphetamine	92 %	95 %	86 %	101 %
	Phentermine	99 %	93 %	90 %	97 %



Extraction of Camphorsulfonic Acid from Serum			
CARTRIDGE	Silia <i>PrepX</i> WAX 3 mL / 60 mg Part Number: SPE-P0020-03BB		
SAMPLE PRETREATMENT	5 mL of serum was mixed with 5 µL of a solution of Camphorsulfonic Acid (0.5 mg/mL) and 5 mL of Phosphoric Acid 4 %		
CONDITIONNING STEP	2 mL of Methanol		
EQUILIBRATION STEP	2 mL of Water		
LOADING STEP	2 mL of sample solution was slowly aspirated through the cartridge		
WASHING STEP	2 mL of 2 % Formic Acid in Water then 2 mL of Methanol, dry the cartridge		
ELUTION STEP	2 mL of 5 % Ammonia in Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS		
RECOVERY	at 0.25 µg/mL: 99 %		

Extra	ction of Alkaloids fro	m Serum
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / Part Number: SPE-P0	3
CONDITIONNING STEP	6 mL of Methanol	
EQUILIBRATION STEP	6 mL of Water	
LOADING STEP	1 mL of serum sample through the cartridge	was slowly aspirated
WASHING STEP	6 mL of Methanol, dry	the cartridge
ELUTION STEP	2 x 3 mL of Acetone	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS	
RECOVERY	at 2 μg/mL	
	Atropine	99 %
	Papaverine	97 %
	Noscapine	95 %
	Strychnine	94 %
	Quinine	60 %

4		Extraction of Caffei	ne, Cotinine & Nicoti	ine from Human L	Irine
CARTRIDGE	Silia <i>PrepX</i> HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
SAMPLE PRETREATMENT	500 μL of urine was	s mixed with 1.5 mL of	Sodium Hydroxide 0.1	М	
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of Water				
LOADING STEP	1 mL of urine sampl	1 mL of urine sample was slowly aspirated through the cartridge			
WASHING STEP	3 mL of Water and dry the cartridge				
ELUTION STEP	3 mL of Methanol				
FURTHER TREATMENT	Evaporation under	Nitrogen, reconstitution	n with Methanol / Water	and quantification b	y LC/MS
RECOVERY	at 100 ng/mL	Silia <i>PrepX</i> HLB	Bond Elut® Plexa	Oasis® HLB	Strata™-X
	Caffeine	97 %	99 %	96 %	97 %
	Cotinine	99 %	100 %	98 %	99 %
	Nicotine	89 %	86 %	90 %	89 %
CONCLUSION	Silia <i>PrepX</i> HLB is	as efficient as competit	tive products to extract (caffeine, cotinine an	d nicotine from urine.



FOOD APPLICATIONS



Extraction of Fungicides in Apple Juice				
CARTRIDGE	Silia <i>PrepX</i> SCX 6 mL / 200 mg Part Number: SPE-P0005-06G			
SAMPLE PRETREATMENT	0.5 mL of Sodium Hyd 5 mL of apple juice	roxide 0.1 N was added to		
CONDITIONNING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of Ammonia 2 %			
LOADING STEP	Treated sample was s cartridge	lowly aspirated through the		
WASHING STEP	3 mL of Ammonia 2 %, 3 mL of 30 % Methanol in Ammonia 5 %, 3 mL of Hydrochloric Acid 0.1 N and 3 mL of Methanol, dry the cartridge			
ELUTION STEP	6 mL of 30 % Methano	ol in Ammonia 5 %		
FURTHER TREATMENT	l '	rogen, reconstitution with quantification by LC/MS		
RECOVERY	at 1 µg/mL			
	Carbendazime	89 %		
	Thiabendazole	92 %		

Extraction of Patulin from Apple Juice			
CARTRIDGE	Silia <i>PrepX</i> HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB		
SAMPLE PRETREATMENT	100 μL of internal standard (250 μg/mL of Patulin-13C (3) in water) and 75 μL of Pectinase Aspergillus Aculeatus were added to 9 mL of apple juice, centrifugate at 3,000 rpm		
CONDITIONNING STEP	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of Water		
LOADING STEP	2 mL of sample supernatant was slowly aspirated through the cartridge		
WASHING STEP	3 mL of 1 % Sodium Bicarbonate and 1 mL of 0.1 % Acetic Acid, dry the cartridge		
ELUTION STEP	2 x 1.5 mL of Ethyl Acetate		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS		
RECOVERY	at 150 ng/kg: 85 %		

	De	ermination of Carbendazim in Orange Juice	
CARTRIDGE	Silia <i>PrepX</i> SCX 3 mL / 60 mg Part Number: SPE-P0005-03B	В	
SAMPLE PRETREATMENT	Centrifugate 5 mL of orange ju Add 2 mL of Acetic Acid 10 % a	ce 5 min at 3,000 rpm. Sample 1 mL of the supernatant. and vortex 1 min	
CONDITIONNING STEP	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of Acetic Acid 10 %		
LOADING STEP	3 mL of the treated sample was	slowly aspirated through the cartridge	
WASHING STEP	2 mL of Acetic Acid 10 % then 2 mL of Methanol		
ELUTION STEP	3 mL of 5 % Ammonium Hydro	kide in Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, re	constitution with Water / Methanol and quantification by LC/MS	
RECOVERY	at 100 ng/mL		
	SiliaPrepX SCX	93 %	
	Bond Elut® Plexa PCX	22 %	
	Oasis® MCX	22 %	
	Strata [™] -X-C	11 %	
CONCLUSION	Silia <i>PrepX</i> SCX performs as w	ell as Waters, Phenomenex & Agilent products to extract carbendazim from orange juice.	





FOOD APPLICATIONS

Enrichment of Streptomycin in Honey			
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G		
SAMPLE PRETREATMENT	Add 2 g of honey to 8 mL of phosphate buffer (pH 2.0), filter, dilute to 16 mL (with the same phosphate buffer) and adjust pH value to 7.5		
CONDITIONNING STEP	5 mL of Methanol		
EQUILIBRATION STEP	3 mL of Water		
LOADING STEP	Treated sample was slowly aspirated through the cartridge		
WASHING STEP	5 mL of Water, dry the cartridge		
ELUTION STEP	5 mL of 3 % Formic Acid in Methanol		
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Water / Acetonitrile and identification by LC/UV		
RECOVERY	at 10 μg/kg: 30 %		

Extra	action of Melamine from Milk
CARTRIDGE	Silia <i>PrepX</i> SCX 6 mL / 200 mg Part Number: SPE-P0005-06G
SAMPLE PRETREATMENT	1 mL of Hydrochloric Acid 1 N was added to 10 mL of milk sample, then mixed with 10 mL of Methylene Chloride. After 15 min centrifugation, remove aqueous layer and extract again organic layer 2 times with 5 mL of Hydrochloric Acid 0.1N. Combine the 3 aqueous fractions.
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of Water
LOADING STEP	Combined aqueous fractions were slowly aspirated through the cartridge
WASHING STEP	6 mL of Hydrochloric Acid 0.1 N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	2x6 mL of 30 % Methanol in Ammonia 5 %
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Water / Methanol and quantification by LC/MS
RECOVERY	at 1 μg/mL: 99 %

	Sulfonamide	s, Tetracyclines &	Pyrimethamine De	etermination in Milk		
CARTRIDGES	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03		lia <i>PrepX</i> DVB 3 mL / art Number: SPE-P00			
SAMPLE PRETREATMENT	Vortex 2 min 600 µL of bovine Adjust pH of the solution at 5.					r (vortex 3 min).
CONDITIONNING STEP	3 mL of Methanol					
EQUILIBRATION STEP	3 mL of Water					
LOADING STEP	1 mL of the treated sample wa	as slowly aspirated th	nrough the cartridge			
WASHING STEP	2 x 3 mL of 10 % Methanol in	Ammonium Acetate	buffer (pH 5.5), dry th	e cartridge		
ELUTION STEP	3 mL of Methanol					
FURTHER TREATMENT	Evaporation under Nitrogen, r	econstitution with Me	ethanol / Water and q	uantification by LC/MS		
FURTHER TREATMENT RECOVERY	Evaporation under Nitrogen, r	reconstitution with Me	ethanol / Water and q	uantification by LC/MS Bond Elut® Plexa	Oasis® HLB	Strata®-X
			,	,		Strata®-X
	at 1,000 pg/mL	Silia <i>PrepX</i> HLB	Silia <i>PrepX</i> DVB	Bond Elut® Plexa	Oasis® HLB	
	at 1,000 pg/mL Sulfathiazol	SiliaPrepX HLB 84 %	SiliaPrepX DVB	Bond Elut® Plexa 85 %	Oasis® HLB	86 %
	at 1,000 pg/mL Sulfathiazol Sulfadiazine	SiliaPrepX HLB 84 % 90 %	SiliaPrepX DVB 83 % 89 %	85 % 88 %	Oasis® HLB 83 % 87 %	86 % 85 %
	at 1,000 pg/mL Sulfathiazol Sulfadiazine Sulfamethoxypyridazine	SiliaPrepX HLB 84 % 90 % 87 %	SiliaPrepX DVB 83 % 89 % 89 %	85 % 88 % 85 %	Oasis® HLB 83 % 87 % 83 %	86 % 85 % 87 %
	at 1,000 pg/mL Sulfathiazol Sulfadiazine Sulfamethoxypyridazine Sulfamethazole	SiliaPrepX HLB 84 % 90 % 87 % 88 %	SiliaPrepX DVB 83 % 89 % 89 % 84 %	85 % 88 % 85 % 87 %	Oasis® HLB 83 % 87 % 83 % 89 %	86 % 85 % 87 % 82 %
	at 1,000 pg/mL Sulfathiazol Sulfadiazine Sulfamethoxypyridazine Sulfamethazole Sulfamethazine	SiliaPrepX HLB 84 % 90 % 87 % 88 % 83 %	SiliaPrepX DVB 83 % 89 % 89 % 84 % 84 %	85 % 88 % 85 % 87 % 86 %	Oasis® HLB 83 % 87 % 83 % 89 % 86 %	86 % 85 % 87 % 82 % 84 %





	Extraction of Marbofloxacin & Sarafloxacin from Salmon	
CARTRIDGE	Silia <i>PrepX</i> SCX 3 mL / 60 mg Part Number: SPE-P0005-03BB	
SAMPLE PRETREATMENT	Add 2 g of salmon and 15 mL of 3 % H ₃ PO ₄ aqueous solution in a 50 mL tube. Shake the tube in a horizontal position for 15 min. Add 5 mL of hexane and vortex for 2 min. Centrifugate at 3,000 rpm for 5 min. Recuperate the aqueous phase from the gelled organic phase by filtration.	
CONDITIONNING STEP	3 mL of Methanol	
EQUILIBRATION STEP	3 mL of Hydrochloric Acid 1 M and 3 mL of water	
LOADING STEP	3 mL of the filtered sample was slowly aspirated through the cartridge	
WASHING STEP	2 mL of Hydrochloric Acid 2 M then 1 mL of Methanol	
ELUTION STEP	3 mL of 10 % Ammonium Hydroxide in Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Water / Methanol and quantification by LC/MS	
RECOVERY	at 50 ppb	
	Marbofloxacin 97 %	
	Sarafloxacin 87 %	

	Extraction of Clenbuterol and Ractopamine from Beef
CARTRIDGE	Silia <i>PrepX</i> WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB
SAMPLE PRETREATMENT	100 μL of internal standard (250 μg/mL of Ractopamine d-6 and 250 μg/mL of Clenbuterol-d9 in Methanol) were added to 1g of chopped beef. Add 5 mL of 0.2 N Sodium Acetate (pH 5.2) and 50 μL of Beta-Glucuronidase / Arylsufatase. Add 2.5 mL of 0.1 M Perchloric Acid, 2 mL of Phosphoric Acid 4 % in Acetonitrile and 5 mL of 0.5 M Glycine (pH 10.5). Adjust to pH 10.5. Add 10 mL of Acetonitrile, 4 g of MgSO ₄ and 1 g of NaCl. Evaporation and reconstitution with 0.1 M Perchloric Acid.
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Water
LOADING STEP	2 mL of treated sample was slowly spirated through the cartridge
WASHING STEP	1.5 mL of Phosphate buffer 25 mM (pH 7.0), then 3 mL of Water and 1 mL of Methanol
ELUTION STEP	3 mL of Formic Acid 2 % in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS
RECOVERY	at 70 ppb
	Clenbuterol 92 %
	Ractopamine 91 %





FOOD APPLICATIONS

	Extraction of Glycoalkaloids from Potatoes
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G
SAMPLE PRETREATMENT	Extract 3 g of potato powder with 20 mL of Water / Acetic Acid / Sodium Metabisulfite (95:5:0.5). Centrifugate for 10 min and filtrer.
CONDITIONNING STEP	5 mL of Acetonitrile
EQUILIBRATION STEP	5 mL of Water / Acetic Acid / Sodium Metabisulfite (95:5:0.5)
LOADING STEP	10 mL of treated sample was slowly aspirated through the cartridge
WASHING STEP	4 mL of 0.5 % Ammonium Hydroxide, then 4 mL of Water and 4 mL of Acetonitrile / Water (15:85). Dry the cartridge.
ELUTION STEP	5 mL of Acetonitrile / Potassium Dihydrogen Phosphate 10 mM (60:40), pH 7.6
FURTHER TREATMENT	Qualitative analysis by TLC

	Acrylamide Determination in Fried Potato Chips
CARTRIDGES	Step 1: SiliaPrepX HLB 6 mL / 200 mg Part Number: SPE-P0002-06G
	Step 2: SiliaPrepX SCX 3 mL / 60 mg Part Number: SPE-P0005-03BB
SAMPLE PRETREATMENT	Extraction 1: vortex for 1 min 1g of potato chips and 8 mL of Sodium Chloride aqueous 4 M. Incubate 30 min at 60°C (vortex 10 sec every 10 min). Centrifugate for 10 min at 4,500 rpm and collect the supernatant.
	Extraction 2: repeat previous 3 steps with same potato chips. Add 1 mL of solution Cirraz 1 (15 g of K_g [Fe(CN)] in 100 mL Water) and 1 mL of solution Cirraz 2 (30 g of $Zn(O_2CCH_3)_2$ in 100 mL Water).
CONDITIONNING STEP (1) [SILIAPREPX HLB]	3 mL of Methanol
EQUILIBRATION STEP (1)	3 mL of Water
LOADING STEP (1)	1.5 mL of the treated sample was slowly aspirated through the cartridge
WASHING STEP (1)	1.5 mL of Water
ELUTION STEP (1)	3 mL of Methanol
CONDITIONNING STEP (2) [SILIAPREPX SCX]	3 mL of Methanol
LOADING STEP (2)	The treated sample eluted from SiliaPrepX HLB was slowly aspirated through the cartridge (collect this fraction)
WASHING STEP (2)	1 mL of Methanol (mix this fraction with the one previously collected)
FURTHER TREATMENT	Evaporation to dryness, reconstitution with Water / Methanol and quantification by LC/MS
RECOVERY	at 100 μg/kg: 88 %





Y.		Determin	ation of Pestici	des in Drinking	Water	
CARTRIDGE	· '	Silia <i>PrepX</i> HLB 6 mL / 200 mg Part Number: SPE-P0002-06G				
CONDITIONNING STEP	6 mL of Methanol					
EQUILIBRATION STEP	6 mL of Water (HPLC	grade)				
LOADING STEP	100 mL of drinking w	ater was slowly as	spirated through tl	ne cartridge		
WASHING STEP	6 mL of Water (HPLC	grade)				
ELUTION STEP	2 x 6 mL of Methanol					
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol and quantification by LC/MS					
RECOVERY	at 1,000 pg/mL	Atrazine	Benalaxyl	Carbendazim	Chloroxuron	lmazalil
	Silia <i>PrepX</i> HLB	75 %	76 %	103 %	91 %	78 %
	Oasis® HLB	66 %	48 %	103 %	99 %	78 %
		Methalaxyl	Myclobutanil	Propoxur	Simazine	Thiambazole
	Silia <i>PrepX</i> HLB	87 %	91 %	70 %	98 %	91 %
	Oasis® HLB	61 %	101 %	42 %	79 %	80 %
CONCLUSION	Silia <i>PrepX</i> HLB comp	pares favorably w	ith Oasis® HLB for	the extraction of	8 pesticides out o	f 10.

Pesticides Determination in Drinking Water				
CARTRIDGE	Silia <i>Prep</i> CleanENVI 6 mL / 1 g Part Number: SPEC-R31930B-06S			
CONDITIONNING STEP	10 mL of Methanol			
EQUILIBRATION STEP	10 mL of Water (HPLC gra	de)		
LOADING STEP	10 mL of drinking water wa through the cartridge	10 mL of drinking water was slowly aspirated through the cartridge		
WASHING STEP	2 x 5 mL of Water (HPLC grade)			
ELUTION STEP	2 x 3 mL of Acetone			
FURTHER TREATMENT	Evaporation under Nitroger Water / Methanol and quar			
RECOVERY	at 50 ng/mL			
	Atrazine	84 %		
	Simazine	95 %		
	Alachlor	68 %		

Pesticides Determination in Water			
CARTRIDGE	Silia <i>PrepX</i> LRV SAX 10 mL / 60 mg Part Number: SPC-P0010-10BB		
CONDITIONNING STEP	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of Water (HPLC grade)		
LOADING STEP	100 mL of sample water was slowly aspirated through the cartridge		
WASHING STEP	3 mL of Water (HPLC grade)		
ELUTION STEP	3 mL of Methanol then 3 mL of Methanol with Formic Acid 2 $\%$		
FURTHER TREATMENT	Quantification by LC/MS/MS		
RECOVERY	at 1,000 pg/mL: > 80 % for 23 pesticides		





Glyphosate & AMPA Determination in Water		
CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB	
SAMPLE PRETREATMENT	Derivatization with FMOC water was added 325 µL 200 µL of EDTA 0.1 M, 4 and 0.6 mL of FMOC-CI Aqueous supernatant wa Acetate. Adjust pH of the adding 100 µL of Formic	of Sodium Borate 50 mM, .5 mL of Acetonitrile 50 mg/mL. Evaporate. as mixed with 2 mL Ethyl a aqeous layer to 3 by
CONDITIONNING STEP	3 mL of Methanol	
EQUILIBRATION STEP	3 mL of Water and 3 mL	of Formic Acid 0.1 %
LOADING STEP	Derivatized sample was the cartridge	slowly aspirated through
WASHING STEP	HING STEP 1 mL of Formic Acid 0.1 % then 2 x 500 μL of Water, dry the cartridge	
ELUTION STEP	TEP 3 mL of Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Water / Acetonitrile and quantification by LC/MS/MS	
RECOVERY	at 5 ng/mL	
	Glyphosate	120 %
	AMPA	106 %

Diquat & Paraquat Determination in Water				
CARTRIDGE	Silia <i>PrepX</i> WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB	3		
CONDITIONNING STEP	3 mL of Methanol			
EQUILIBRATION STEP	3 mL of Water			
LOADING STEP	100 mL of sample water was slowly aspirated through the cartridge			
WASHING STEP	3 mL of Water then 3 mL of Methanol			
ELUTION STEP	2 x 3 mL of Acetonitrile / Isopropanol / Formic Acid (40:40:20)			
FURTHER TREATMENT	Quantification by LC/MS/MS			
RECOVERY	at 10 ppb			
	Diquat 9	0 %		
	Paraquat 9	0 %		

		Determinat	tion of Pesticides in W	ater (<i>by GC/ECI</i>	D)
CARTRIDGE	Silia <i>PrepX</i> HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
CONDITIONNING STEP	3 mL of 30 % Acetone in 1	Toluene then	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of distilled Water				
LOADING STEP	100 mL of sample water wa	100 mL of sample water was slowly aspirated through the cartridge			
WASHING STEP	3 mL of distilled Water, dry the cartridge				
ELUTION STEP	$500~\mu L$ of Acetone, then 2 mL of 30 % Acetone in Toluene and 2.5 mL of 30 % Acetone in Toluene				
FURTHER TREATMENT	Qualification by GC/ECD				
RECOVERY	Trifluralin	90 %	Endrin	95 %	
	Lindane	88 %	4,4'-DDT	75 %	
	Aldrin	78 %	Diclofop-methyl	90 %	
	Heptachlor Epoxide	88 %	Methoxychlor	86 %	
	Dieldrin	90 %	Chlordane	79 %	





Extraction of Desphenyl Chloridazon from Water		
CARTRIDGE	Silia <i>PrepX</i> SAX 3 mL / 60 mg Part Number: SPE-P0010-03BB	
SAMPLE PRETREATMENT	100 μL of Ammonium Hydroxide 5 % was added to 1 mL of water sample	
CONDITIONNING STEP	1 mL of Methanol	
EQUILIBRATION STEP	1 mL of Ammonium Hydroxide 5 %	
LOADING STEP	Treated sample was slowly aspirated through the cartridge	
WASHING STEP	1 mL of Ammonium Hydroxide 5 % then 1 mL of Methanol	
ELUTION STEP	2 x 1 mL of 5 % Formic Acid in Ethyl Acetate	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS	
RECOVERY	at 10 μg/mL: 83 %	

Quantification of Acidic Herbicides				
CARTRIDGE	Silia <i>PrepX</i> SAX 6 mL / 200 mg Part Number: SPE-P0010-06G			
SAMPLE PRETREATMENT	pH of sample was adjusted to basic value with Sodium Hydroxide 1 N			
CONDITIONNING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of Water			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	3 mL of Sodium Acetate then 3 mL of Methanol, dry the cartridge			
ELUTION STEP	2 x 3 mL of Formic Acid 10 % in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS			
RECOVERY	at 1 μg/mL			
	Bentazon	79 %		
	Dicamba	87 %		
	2,4-Dichlorophenoxy Acetic Acid	82 %		
	4-Chloro-2-methylphenoxy Acetic Acid	76 %		

	Isothi	azolinone Biocides in an Aqueous Sample	
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-06G		
SAMPLE PRETREATMENT	5 mL of Isothiazolinones standard solution (1 μ g/mL) are diluted in 50 mL Water and 500 μ L Formic Acid. The solution is filled up to 100 mL.		
CONDITIONNING STEP	6 mL of Methanol		
EQUILIBRATION STEP	6 mL of 0.1 % Formic Acid in Wat	er	
LOADING STEP	5 mL of sample was slowly aspirate	ed through the cartridge	
WASHING STEP	6 mL 0.1 % Formic Acid in Water, dry the cartridge		
ELUTION STEP	3 mL of Methanol then 6 mL of Ac	cetonitrile	
FURTHER TREATMENT	Evaporation under Nitrogen, reco	onstitution with Methanol / Water and quantification by LC/MS	
RECOVERY	at 50 ng/L		
	Methylisothiazolinone	93 %	
	Chloromethylisothiazolinone	96 %	
	Benzisothiazolinone	85 %	
	Butylbenzisothiazolinone	88 %	
	Octylisothiazolinone	90 %	
	Dichloroctylisothiazolinone	83 %	





	Extraction of Bisphenol A, Triclosan & Ethynyl Estradiol from Water
CARTRIDGE	Silia <i>PrepX</i> HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB
SAMPLE PRETREATMENT	To 25 mL of sample water was added 250 μL of internal standard (1 ppb of 17α-Ethynyl Estradiol d-6, 1 ppb of Bisphenol A d-16 and 0.4 ppb of Triclosan d-3 in Methanol)
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Water and 1 mL of Acetic Acid 100 mM
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	3 mL of Water, 1 mL of Acetic Acid 100 mM and 2 mL of 20 % Methanol in Water, dry the cartridge
ELUTION STEP	2 x 3 mL of Dichloromethane / Acetone (50:50)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Sodium Carbonate in Water, derivatization with Dansyl Chloride and quantification by LC/MS/MS
RECOVERY	17α-Ethynyl Estradiol 93 %
	Bisphenol A 115 %
	Triclosan 98 %

	Analysis of Bisphenol A in Bottled Water
CARTRIDGE	Silia <i>Prep</i> C18 Plus 6 mL / 200 mg (<i>glass</i>) Part Number: SPE-R00830B-06G
SAMPLE PRETREATMENT	100 μL of internal standard (<i>Bisphenol A-d16 in methanol, 1 μg/mL</i>) was added to 50 mL of bottled water
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Water (HPLC grade) and 1 mL of Acetic Acid 100 mM
LOADING STEP	The whole sample was aspirated through the cartridge using SiliCycle MiniBlock equipment (2 drops / second)
WASHING STEP	5 mL of Water (<i>HPLC grad</i> e), dry the cartridge
ELUTION STEP	3 mL of Methanol
FURTHER TREATMENT	Evaporation to dryness, derivatization using Dansyl Chloride, liquid-liquid extraction, evaporation, reconstitution with Methanol and quantification by LC/MS/MS
RECOVERY	at 3,000 pg/mL: 97 %





		Pharmac	eutical Drugs Determi	nation in Water		
CARTRIDGE	Silia <i>PrepX</i> HLB (200 mg) + SAX (60 mg) / 10 mL Part Number: SPC-P0210-10i					
SAMPLE PRETREATMENT			ith 5 mL of Sodium Acetat olution (<i>NH₄Cl 0.5 M and l</i>			
CONDITIONNING STEP	6 mL of Methanol					
EQUILIBRATION STEP	6 mL of Water and 6 mL of	of buffer pH 9	0.5			
LOADING STEP	Treated sample was slow	ly aspirated t	through the cartridge			
WASHING STEP	3 mL of buffer pH 9.5 and	3 mL of buffer pH 9.5 and 3 mL of Water, dry the cartridge				
ELUTION STEP	2 mL of Methanol and 2 r	2 mL of Methanol and 2 mL of Formic Acid 2 % in Methanol				
FURTHER TREATMENT	Evaporation under Nitrog	Evaporation under Nitrogen, reconstitution with Water / Acetonitrile and quantification by LC/MS/MS				
RECOVERY	at 100 ppt					
	Trimethroprim	105 %	Caffeine C13	96 %		
	Sulphamethoxazole	100 %	Acetaminophen	93 %		
	Naproxen	Naproxen 100 % Norfloxacin 70 %				
	Ibuprofen	85 %	Maprotiline	79 %		
	Carbamazepine	102 %				

Determination of Tricyclic Antidepressants in Water			
CARTRIDGE	Silia <i>PrepX</i> DVB 3 mL / 60 mg Part Number: SPE-P0001-03BB		
CONDITIONNING STEP	1 mL of Methanol		
EQUILIBRATION STEP	1 mL of Water		
LOADING STEP	1 mL of sample was slowly aspirated through the cartridge		
WASHING STEP	1 mL of Water		
ELUTION STEP	1 mL of Acetonitrile		
FURTHER TREATMENT	Quantification by LC/MS		
RECOVERY	at 1 μg/mL		
	Protriptyline	93 %	
	Nortriptyline	90 %	





Y.		Determin	ation of Explosives in Well W	ater
CARTRIDGE	Silia <i>PrepX</i> DVB 6 mL / 200 mg Part Number: SPE-P0001-060			
CONDITIONNING STEP	6 mL of Methanol, 6 mL of Ace	etonitrile		
EQUILIBRATION STEP	10 mL of Water			
LOADING STEP	1 L of well water (with 5 g of S	odium Chlo	oride) was slowly aspirated throug	h the cartridge
WASHING STEP	10 mL of Water, DO NOT dry t	the cartridg	e	
ELUTION STEP	6 of mL Methanol / Acetonitrile	e (50:50)		
FURTHER TREATMENT	Evaporation under Nitrogen, r	econstitutio	on with Methanol / Water and quar	ntification by LC/MS
RECOVERY	at 1 µg/L			
	Hexanitrodiphenylamine	96 %	4-Amino-2,6-dinitrotoluene	95 %
	Diphenylamine	100 %	2-Amino-4,6-dinitrotoluene	94 %
	Pentaerythritol Tetranitrate	108 %	2,4,6-Trinitrotoluene	92 %
	3-Nitrotoluene	78 %	Nitroglycerine	88 %
	4-Nitrotoluene	81 %	1,3-Dinitrobenzene	86 %
	2-Nitrotoluene	67 %	1,3,5-Trinitrobenzene	96 %
	2,6-Dinitrotoluene	94 %	Ethylene Glycol Dinitrate	95 %
	2,4-dinitrotoluene	85 %	Picric Acid	92 %
	Octogen	94 %	Diethylene Glycol Dinitrate	74 %

Determination of Surfactants in Water				
CARTRIDGE	Silia <i>PrepX</i> WAX 3 mL / 60 mg Part Number: SPE-P0020-03BB			
CONDITIONNING STEP	2 mL of 5 % Ammonia in Methanol then 2 mL of Methanol			
EQUILIBRATION STEP	2 mL of Water			
LOADING STEP	500 mL of water sample was slowly aspirate through the cartridge	d		
WASHING STEP	2 mL of Water, then 2 mL of Acetone / Aceton Formic Acid (50:50:1) and 2 mL of Methanol	itrile /		
ELUTION STEP	2 mL of 5 % Ammonia in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution w Methanol / Water and quantification by LC/N			
RECOVERY	at 20 µg/L			
	Perfluorooctane Sulfonate Potassium Salt	81 %		
	Perfluoropentanoic Acid	94 %		
	Perfluorohexanoic Acid	94 %		
	Perfluorooctanoic Acid	95 %		
	Perfluoropropionic Acid	103 %		
	Perfluorododecanoic Acid	82 %		

Quantification of Phenolic Acids					
CARTRIDGE	SiliaPrepX SAX 6 mL / 200 mg Part Number: SPE-P0010-06G				
SAMPLE PRETREATMENT	pH of sample was adjusted to basic value with Sodium Hydroxide 1 N				
CONDITIONNING STEP	6 mL of Methanol				
EQUILIBRATION STEP	6 mL of water				
LOADING STEP	Treated sample was slowly aspirated through the cartridge				
WASHING STEP	3 mL of water, then 3 mL of Sodium Hydroxide 0.1 N and 3 mL of Methanol, dry the cartridge				
ELUTION STEP	2 x 3 mL of Formic Acid 5 % in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen Methanol / water and quant	,			
RECOVERY	at 1 μg/mL				
	Syringic Acid	70 %			
	Vanillic Acid 86 %				
	p-Hydroxybenzoic Acid 97 %				





	Extraction	on of Amines	from an Aqueous Sample				
CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G						
SAMPLE PRETREATMENT	200 μL of Phosphoric Acid 2 % was a	dded to 1 mL	of aqueous sample				
CONDITIONNING STEP	6 mL of Methanol						
EQUILIBRATION STEP	6 mL of Water						
LOADING STEP	Treated sample was slowly aspirated	through the ca	artridge				
WASHING STEP	6 mL of Hydrochloric Acid 0.1 N then	6 mL of Hydrochloric Acid 0.1 N then 6 mL of Methanol, dry the cartridge					
ELUTION STEP	2 x 3 mL of 10 % Ammonia in Methanol						
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / Water and quantification by LC/MS						
RECOVERY	at 100 ppm						
	2-Naphthylamine	65 %	4,4'-Methylene-bis-(2-chloro-aniline)	75 %			
	Benzidine	104 %	4,4′-Oxydianiline	104 %			
	5-Nitro-o-toluidine	80 %	4,4´-Methylenedianiline	109 %			
	Xenylamine	89 %	4,4´-Thiodianiline	100 %			
	o-Aminoazotoluene	o-Aminoazotoluene 89 % 4,4´-Methylendi-o-toluidine 110 %					
	4-Aminoazobenzene	99 %	3,3-Dichlorobenzidine	110 %			

Y.		Extraction of PAHs f	rom Drinking Water			
CARTRIDGE	SiliaPrep PAH 6 mL / 1.5 g Part Number: SP2-R0610030B-06T					
CONDITIONNING STEP	5 mL of 2-Propanol	5 mL of 2-Propanol				
EQUILIBRATION STEP	5 mL of water / 2-Propanol (9	92:8)				
LOADING STEP	500 mL of drinking water was	s slowly aspirated through	the cartridge			
WASHING STEP	3 mL of Dichloromethane HP	LC grade, soak the sorber	nt for 10 minutes before elutir	ng. Repeat a second time.		
ELUTION STEP	2 mL of Dichloromethane HP	LC grade, soak the sorber	nt for 10 minutes before elutir	ng. Combine the 3 eluates.		
FURTHER TREATMENT	Evaporation under Nitrogen,	reconstitution with Acetone	e / Water and qualification by	HPLC (Fluorescence)		
RECOVERY		Benzo[b]fluoranthene	Benzo[k]fluoranthene	Benzo[a]pyrene		
	Silia <i>Prep</i> PAH	118 %	99 %	94 %		
	BAKERBOND PAH Aqua	117 %	102 %	100 %		
		Benzo[ghi]perylene	Indeno[1,2,3-cd]pyrene			
	Silia <i>Prep</i> PAH	117 %	126 %			
	BAKERBOND PAH Aqua	115 %	114 %			
CONCLUSION	Silia <i>Prep</i> PAH performs as w	rell as BAKERBOND PAH	Aqua for the extraction of PA	Hs from water.		





	Analysis of Pesticides in Oats, after a Fatty Acids Cleanup
CARTRIDGE	Silia <i>Prep</i> Diamine 6 mL / 500 mg Part Number: SPE-R49030B-06P
SAMPLE PRETREATMENT	10 g of oat was added to 100 mL of Water and 200 mL of Acetone. 35 g of NaCl and 100 mL of 50 % Ethylacetate in Cyclohexane were added for liquid-liquid extraction. The organic layer (200 mL) was dried with NaSO ₄ , filtered, evaporated and reconstituted with 10 mL of 50 % Ethylacetate in Cyclohexane.
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Acetone and 3 mL of 50 % Ethyl Acetate in Cyclohexane
LOADING STEP	1 mL of treated sample was slowly aspirated through the cartridge (collect the eluted solvent)
WASHING STEP	6 mL of Hydrochloric Acid 0.1 N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	15 mL of 50 % Ethyl Acetate in Cyclohexane (mix with the fraction previously collected)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile, derivatization using HMDS and TFA, and quantification by GC/MS
RECOVERY	> 80 % for 84 pesticides
	< 1 % for fatty acids

Source: Steinbach P. et al., J. Chromatogr. A, 2014, 1323, 28 - 38

Triacylglycerols Profiling of Marine Microalgae		
CARTRIDGE	Silia <i>Prep</i> Silica 3 mL / 500 mg Part Number: SPE-R10030B-03P	
SAMPLE PRETREATMENT	Algal extracts were extracted with Hexane, washed with water and evaporated	
CONDITIONNING STEP	3 mL of Hexane	
EQUILIBRATION STEP	3 mL of distilled Water	
LOADING STEP	50 mg of lipid sample in 300 μL of Hexane was slowly aspirated through the cartridge	
ELUTION STEP	Elution 1 (for triacylglycerols): Hexane / Diethyl Ether / Acetic Acid (80:20:1) Elution 2 (for polar lipids and chlorophyl): Acetone	
FURTHER TREATMENT	Evaporation, reconstitution with Hexane and quantification by LC/MS/MS	

Source: Franz A. et al., Journal of Lipid Research, 2011, 52, 2101 - 2108

Extraction of Allantoin from a Cosmetic Product			
CARTRIDGE	Silia <i>PrepX</i> SAX 3 mL / 60 mg Part Number: SPE-P0010-03BB		
SAMPLE PRETREATMENT	1 g of cosmetic was diluted in 100 mL of Water, pH was adjusted to 10°C with Ammonium Hydroxide 5 %		
CONDITIONNING STEP	3 mL of Methanol		
EQUILIBRATION STEP	3 mL of Ammonium Hydroxide 5 %		
LOADING STEP	1 mL of treated sample was slowly aspirated through the cartridge		
WASHING STEP	3 mL of Ammonium Hydroxide 5 % then 3 mL of Methanol		
ELUTION STEP	2 x 1 mL of Hydrochloric Acid 0.6 %		
FURTHER TREATMENT	Add Acetonitrile / Ammonium Chloride 30 mM, and qualification by HPLC		



Share your complete Silia*Prep I* Silia*PrepX* application with SiliCycle... ... and get any standard* SPE box for free!

Send us:

- The detailed method (cartridge used, conditionning, equilibration, loading, washing & elution steps)
- · Analytes & Matrix
- Pre- and Post-treatments
- Recoveries
- HPLC or GC analysis

Contact us: sampleprep@silicycle.com



^{*} box of 200 cartridges are not available for this promotion

Suggested Protocols for Various Analytes & Matrices

	Chart of Suggested Protocols							
Application	Matrix	Analytes	SPE Cartridge	Sample Pretreatment before Loading				
Aflatoxins	Peanut butter	Aflotoxin M1, Aflotoxin G2, Aflotoxin G1, Aflotoxin B2, Aflotoxin B1	SiliaPrep Florisil PR 500 mg / 3 mL SPE-AUT-0015-03P and SiliaPrep Silica 200 mg / 3 mL SPE-R10030B-03G	Add 40 mL of Methanol / Water (80:20) and 0.2 g of Sodium Chloride to 5 g of peanut butter. Stir for 2 hours, filter on paper and rinse with 15 mL of Methanol. Dry the combined extracts with Magnesium Sulfate and evaporate under Nitrogen. Reconstitute with 500 µL of Methanol / Water (80:20).				
Aminoglycoside Antibiotics	Milk, meat (beef, chicken, pork) and eggs	Spectinomycin, Apramycin, Dihydrostreptomycin, Gentamicin, Hygromycin B, Kanamycin, Neomycin B, Streptomycin, Amikacin, Netilmicin, Paromomycin, Sisomycin, Tobramycin	Silia <i>PrepX</i> WCX 500 mg / 6 mL SPE-P0015-06P	Weigh 5 g of mixed eggs or frozen and milled meat samples, or 10 mL of bovine milk. Add 20 mL of the extraction solution: NH4OAC (10 mM), EDTA (0.4 mM), NaCl (0.5 %) and Trichloroacetic Acid (2 %) in Water. Vortex and centrifugate at 4,000 rpm for 10 minutes. Transfer the supernatant to a clean tube. Repeat extraction and combine supernatants. Adjust pH to 6.5.				
Barbiturates	Biological fluids (blood, plasma, serum, urine or tissue)	Phenobarbital, Butalbital, Amobarbital, Pentobarbital, Secobarbital, Butabarbital, Hexobarbital	Silia <i>Prep</i> CleanDRUG 200 mg / 6 mL SPEC-R651230B-06G	Add 3 mL of Phosphate buffer (100 mM, pH 6.0) to 2 mL of blood / plasma / serum / urine (or 1 g of tissue homogenate). Mix / vortex. Adjust pH to 6.0 ± 0.5 with Sodium Phosphate (100 mM). Centrifugate at 2,000 rpm for 10 minutes and discard cellular debris.				
Benzodiazepines	Urine	Nordiazepam, 7-Aminoclonazepam, Desalkylflurazepam, Temazepam, Alprazolam, Clonazepam, Midazolam, Flurazepam, 7-Aminoflunitrazepam, Diazepam, Oxazepam, Chlordiazepoxide, Flunitrazepam, Lorazepam	SiliaPrepX SCX 30 mg / 3 mL SPE-P0005-03AA	Dilute 1 mL of urine with 2 mL of 2 % Formic Acid.				
Biomarkers of Alcohol Consumption	Urine	Ethyl Glucuronide (<i>EtG</i>), Ethyl Sulfate (<i>EtS</i>)	SiliaPrep Amine (WAX) 100 mg / 1 mL SPE-R52030B-01C	Add 50 µL of HCl (6 M) and 1 mL of Acetonitrile to 100 µL of urine, centrifugate.				
Buprenorphine & Norbuprenorphine	Whole blood	Buprenorphine, Norbuprenorphine	SiliaPrepX SCX 30 mg / 3 mL SPE-P0005-03AA	Dilute 0.5 mL of blood with 3 mL of Phosphate buffer (0.1 M, pH 6.0). Vortex, centrifugate and discard cellular debris.				
Caffeine, Theophylline & Theobromine	Biological fluids (blood, plasma, serum, urine or tissue)	Caffeine, Theophylline, Theobromine	Silia <i>Prep</i> CleanDRUG 200 mg / 6 mL SPEC-R651230B-06G	Add 3 mL of Acetic Acid (100 mM) to 1 mL of blood / serum / plasma / urine. Vortex, centrifugate.				
Cocaine and Metabolites	Biological fluids (whole blood, urine or plasma)	Ecgonine Methyl Ester, Anhydroecgonine Methyl Ester, Benzoylecgonine, Norcocaine, Cocaine, Cocaethylene	Silia <i>PrepX</i> SCX 1 mL / 30 mg SPE-P0005-01AA	Plasma and Urine: dilute 100 μL of sample with 300 μL of Ammonium Acetate buffer (50 mM, pH 6.0). Whole blood: lyse red blood cells by 10 minutes sonication in buffer, centrifugate at 11,000 rpm for 10 minutes, discard cellular debris.				
Explosives	Surface water	Nitrobenzene, 2-Nitrotoluene, 3-Nitrotoluene, 4-Nitrotoluene, 1,3-Dinitrobenzene, 2,6-Dinitrotoluene, 2,4-Dinitrotoluene, 1,3,5-Trinitrobenzene, 2,4,6-Trinitrotoluene, RDX, 4-Amino-2,6-Dinitrotoluene, 3,5-Dinitroaniline, 2-Amino-4,6-Dinitrotoluene, Tetryl	SiliaPrepX DVB 500 mg / 6 mL SPE-P0001-06P	No sample pretreatment required.				

Note: convenient starting points for method development, further optimization may be required



Chart of Suggested Protocols							
Con	ditionning	Equilibration	Loading	Washing	Elution	Further Treatment	
when	er recoveries n no itionning	Better recoveries when no equilibration	Pretreated sample	2 x 3 mL of Methanol / Water (80:20) 2 x 3 mL of Methanol	2 x 3 mL of Acetone / Water / Formic Acid (96:3.5:0.5)	Combine elution fractions, evaporate under Nitrogen and reconstitute with 2 mL of Hexane / Chloroform (50:50). Condition the SiliaPrep Silica cartridge with 2 x 3 mL of Hexane and load the sample extracted from SiliaPrep Florisil PR cartridge. Wash with 3 x 2 mL of Methanol / Chloroform (50:50). Combine loading and washing fractions, evaporate under Nitrogen, reconstitute with 0.1 % Formic Acid and Ammonium Acetate (5 mM) in Water / Methanol (50:50) and quantify by LC/MS.	
10 m Meth	-	10 mL of Water	Pretreated sample	10 mL of Ammonium Acetate buffer (<i>50 mM, pH 7.0</i>) Dry the cartridge	6 mL of 2 % Formic Acid in Water	Dilute with 4 mL of Water, filter on 0.2 µm and quantification by LC/MS/MS.	
3 mL	of Methanol	3 mL of Phosphate buffer (100 mM, pH 6.0)	Pretreated sample	3 mL of Water 1 mL of Acetic Acid (100 mM) Dry the cartridge 2 mL of Hexane	3 mL of Ethyl Acetate / Hexane (50:50)	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water and quantification by LC/MS/MS. Or dissolution in Ethyl Acetate and quantification by GC/MS.	
1 mL	of Methanol	1 mL of Water	Pretreated sample	2 mL of 2 % Formic Acid 2 mL of Methanol / Water (50:50)	1 mL of 5 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Methanol and quantification by LC/MS/MS.	
2 mL	of Methanol	2 mL of Water 2 mL of 0.2 % Acetic Acid in Acetonitrile	Pretreated sample	1 mL of Hexane Dry the cartridge	2 x 750 µL of 2 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with Water and quantification by LC/MS/MS.	
3 mL	of Methanol	3 mL of water	Pretreated sample	2 x 2 mL of 2 % Formic Acid 3 mL of Methanol / 2 % Formic Acid (70:30) Dry the cartridge	2 x 1.5 mL of 5 % Ammonium Hydroxide in Ethyl Acetate / Isopropanol (80:20)	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Methanol / Water (15:85) and quantification by LC/MS/MS.	
	of Methanol of Water	1 mL of Acetic Acid (100 mM)	Pretreated sample	3 mL of Water 3 mL of Acetic Acid (100 mM) Dry the cartridge	3 mL of Ethyl Acetate / Methanol (90:10)	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water and identification by HPLC.	
1 mL	of Methanol	1 mL of Ammonium Acetate (50 mM, pH 6.0)	Pretreated sample	1 mL of Ammonium Acetate buffer (50 mM, pH 6.0) 1 mL of 2 % Formic Acid in Water 1 mL of Methanol	2 x 0.5 mL of 5 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with Water / Methanol (80:20) and quantification by LC/MS.	
6 mL Aceta	of Ethyl ate	6 mL of Methanol 6 mL of Water	1 L of water sample Dry the cartridge	No washing required	6 mL of Ethyl Acetate	Evaporation down to 1 mL under Nitrogen and identification by HPLC.	



Suggested Protocols for Various Analytes & Matrices (continued)

Chart of Suggested Protocols									
Application	Matrix	Analytes	SPE Cartridge	Sample Pretreatment before Loading					
FAMEs	Olive oil	C11:0, C16:0, C16:1 cis 9, C18:0, C18:1 cis 9, C18:2 cis 9,12, C18:3 cis 9,12,15, C20:0, C20:1 cis 11	Silia <i>Prep</i> Silica 1 g / 6 mL SPE-R10030B-06S	Dilute 0.12 g of oil in 0.5 mL of Hexane.					
Fermentable Sugars	Beer	Fructose, Glucose, Maltose, Maltotriose, Maltotetraose, Maltopentaose, Maltohexaose, Maltoheptaose	SiliaPrep C18 Plus 500 mg / 6 mL SPE-R00830B-06P	Remove carbon dioxide by shaking the beer.					
Fluoroquinolones	Milk	Norfloxacin, Ofloxacin, Ciprofloxacin, Pefloxacin, Lomefloxacin, Danofloxacin, Enrofloxacin, Sarafloxacin, Difloxacin, Oxolinic Acid, Flumequine	SiliaPrepX DVB 60 mg / 3 mL SPE-P0001-03BB	Add 15 mL of Trichloroacetic Acid / Acetonitrile (10:90) to 1 g of milk. Mix / vortex. Centrifugate at 5,000 rpm at 4°C for 10 minutes.					
Halocetic Acids	Surface water	Chloroacetic Acid, Bromoacetic Acid, Dichloroacetic Acid, Dalapon, Trichloroacetic Acid, Bromochloroacetic Acid, Bromodichloroacetic Acid, Dibromoacetic Acid, Chlorodibromoacetic Acid, Tribromoacetic Acid	SiliaPrep SAX 500 mg / 6 mL SPE-R66530B-06P	Add 0.5 mL of aqueous Ammonium Chloride to 50 mL of water sample. Adjust pH to 4.5 - 5.5 with Sulfuric Acid.					
Hormones	Serum	Aldosterone, Cortisol, Cortisone, Corticosterone, 11-Deoxycortisol, β -Estradiol, Testosterone, 11-Deoxycorticosterone, Androstenedione, Estrone, 17 α OH Progesterone, DHT, Progesterone	SiliaPrepX DVB 30 mg / 1 mL SPE-P0001-01AA	Add 400 µL of a 0.5 % Formic Acid solution to 100 µL of serum. Vortex, centrifugate, discard cellular debris.					
Lipids	Tissue	Fatty Acids, Phospholipids, Cholesteryl Ester, Triglycerides, Cholesterol, Diglycerides, Monoglycerides	Silia <i>Prep</i> Amine (<i>WAX</i>) 500 mg / 3 mL SPE-R52030B-03P	Evaporate extract to dryness under Nitrogen and dissolve in 0.5 mL of Chloroform.					
Melamine & Analogues	Powdered infant milk	Melamine, Ammeline, Ammelide	SiliaPrepX SCX 60 mg / 3 mL SPE-P0005-03BB	Add 20 mL of Water to 1 g of powdered infant milk, vortex. Take 1 mL and add 2 mL of HCl (0.1 N).					
Melamine & Analogues	Powdered infant milk	Cyanuric Acid	SiliaPrepX SAX 60 mg / 3 mL SPE-P0010-03BB	Add 20 mL of Water to 1 g of powdered infant milk, vortex. Take 1 mL and add 2 mL of Sodium Hydroxide (0.1 N).					
Metal Ions	Surface water	Cu(II), Hg(II), Se(IV), Zn(II)	Silia <i>Prep</i> Amine (<i>WAX</i>) 50 mg / 1 mL SPE-R52030B-01B	No sample pretreatment required.					
MMA & Succinic Acid	Plasma	Methylmalonic Acid (<i>MMA</i>), Succinic Acid	SiliaPrepX WAX 30 mg / 1 mL SPE-P0020-01AA	Dilute 100 µL of plasma with 0.5 mL of 0.1 % Acetic Acid.					
Opiates	Urine	Morphine, Codeine, Hydromorphone, Norcodeine, Hydrocodone, Oxycodone, Oxymorphone	Silia <i>PrepX</i> SCX 30 mg / 3 mL SPE-P0005-03AA	Add 125 µL of concentrated HCl to 0.5 mL of urine. Heat at 95°C for 1.5 hour. Cool, add 2 mL of Sodium Acetate buffer (0.1 M, pH 4.5). Neutralize with 250 µL of Potassium Hydroxide (7 N), vortex. pH should be inferior to 6.0. Centrifugate at 6,000 rpm for 20 minutes.					

Note: convenient starting points for method development, further optimization may be required



Chart of Suggested Protocols						
Conditionning	Equilibration	Loading	Washing	Elution	Further Treatment	
No conditioning required	6 mL of Hexane	Pretreated sample	No washing required	10 mL of Hexane / Diethyl Ether (87:13)	Evaporation under Nitrogen, reconstitution with Hexane, esterification and identification by GC/FID.	
6 mL of Methanol	6 mL of Water	0.5 mL of beer	No washing required	6 mL of Water	Identification by HPLC.	
3 mL of Methanol	3 mL of Water	Pretreated sample	3 mL of Methanol / Water (10:90) Dry the cartridge	3 mL of Methanol	Evaporation under Nitrogen, reconstitution with mobile phase and identification by HPLC.	
10 mL of Methanol	10 mL of Water	Pretreated sample	10 mL of Methanol	3 mL of 10 % Sulfuric Acid / Methanol	Esterification and identification by GC/µECD.	
1 mL of 0.5 % Formic Acid in Methanol	1 mL of 0.5 % Formic Acid in Water	Pretreated sample	1 mL of Methanol / Water (30:70) Dry the cartridge	2 x 0.25 mL of Methanol	Evaporation under Nitrogen, reconstitution with Methanol / Water (50:50) and quantification by LC/MS/MS.	
No conditionning required	2 x 2 mL of Hexane	Pretreated sample	 Elution of Neutral Lipids: 4 mL of Chloroform / 2-Propanol (2:1). Elution of Fatty Acids: 4 mL of 2 % Acetic Acid in Diethyl Ether. Elution of Phospholipids: 4 mL of Methanol. With Neutral Lipids fraction: Evaporation under Nitrogen, reconstitution with 0.2 mL of Hexane, load of a second SPE cartridge (same equilibration step). Elution of Cholesteryl Ester: 4 mL of Hexane. Connect a third SPE cartridge below the second one (same equilibration step). Elution of Triglycerides: 6 mL of Diethyl Ether / Methylene Chloride / Hexane (1:10:89), through both SPE cartridges. Elution of Cholesterol: 12 mL of 5 % Ethyl Acetate in Hexane, through both SPE cartridges. On the second SPE cartridge only now. Elution of Diglycerides: 4 mL of 15 % Ethyl Acetate in Hexane. Elution of Monoglycerides: 4 mL of Chloroform / Methanol (2:1). 			
3 mL of Methanol	3 mL of Water	Pretreated sample	3 mL of 2 % Formic Acid in Water 3 mL of Methanol / Acetonitrile (50:50)	3 mL of 5 % Ammonium Hydroxide in Methanol / Acetonitrile (50:50)	Evaporation under Nitrogen, reconstitution with Methanol / Acetonitrile (50:50) and quantification by LC/MS.	
3 mL of Methanol	3 mL of Water	Pretreated sample	3 mL of Water 3 mL of Methanol	3 mL of 5 % Acetic Acid in Methanol	Evaporation under Nitrogen, reconstitution with Methanol / Acetonitrile (50:50) and quantification by LC/MS.	
3 mL of Methanol	3 mL of Water	50 mL water sample	3 mL of Water Dry the cartridge	3 mL of a Nitric Acid solution (100 mM)	Quantification by ICP-AES.	
1 mL of Methanol	1 mL of 0.1 % Acetic Acid in Water	Pretreated sample	0.5 mL of Methanol / Water (50:50) Dry the cartridge	2 x 0.6 mL of 2 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water and quantification by LC/MS/MS.	
0.5 mL of Methanol	1 mL of Water	Pretreated sample	1 mL of 2 % Formic Acid 1 mL of Methanol Dry the cartridge	2 mL of Ammonium Hydroxide / Methanol (20:100)	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Methanol / Water (5:95) and quantification by LC/MS/MS.	



Suggested Protocols for Various Analytes & Matrices (continued)

	Chart of Suggested Protocols					
Application	Matrix	Analytes	SPE Cartridge	Sample Pretreatment before Loading		
Organophosphates	Surface water	Tris (1-Chloro-2-Propyl)-Phosphate (<i>TCPP</i>), Tris (2-Chloroethyl)-Phosphate (<i>TCEP</i>), Tris (1,3-Dichloro-2-Propyl)-Phosphate (<i>TDCP</i>), Tri-n-Butylphosphate (<i>TnBP</i>), Tri-Isobutylphosphate (<i>TiBP</i>), Tris(2-Butoxyethyl)-Phosphate (<i>TBEP</i>)	SiliaPrepX HLB 100 mg / 1 mL SPE-P0002-01C	No sample pretreatment required.		
PAHs	Chocolate	Naphthalene, Acenaphthylene, Acenaphthene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benzo[c]fluorene, Benzo[a]anthracene, Chrysene, Cyclopental[cd]pyrene, Triphenylene, 5-methyl-Chrysene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Dibenz[a,h]anthracene, Indeno[1,2,3-cd]pyrene, Benzo[ghi]perylene, Dibenzo[a,l]pyrene, Dibenzo[a,e]pyrene, Dibenzo[a,i]pyrene, Dibenzo[a,h]pyrene	Silia <i>Prep</i> Silica 1 g / 6 mL SPE-R10030B-06S	Grind 1 g of chocolate with 10 mL of Methanol to extract fats and coca butter. Evaporate Methanol under Nitrogen. Add 10 mL of Water then 5 mL of n-Pentane and proceed to liquid-liquid extraction. Repeat the extraction with another 5 mL of n-Pentane. Combine the two portions of 5 mL of n-Pentane, evaporate to 2 mL under Nitrogen.		
Parabens	Cosmetics (toothpaste)	Methyl Paraben, Propyl Paraben	Silia <i>Prep</i> C8 500 mg / 3 mL SPE-R31030B-03P	Add 10 mL of Methanol to 1 g of toothpaste, vortex, centrifugate. Take 100 μL of mixture and dilute to 2 mL with Methanol.		
РСВ	Surface water	2-Chlorobiphenyl, 4-Chlorobiphenyl, 2,4'-Dichlorobiphenyl, 2,2',5-Trichlorobiphenyl, 2,4,4'-Trichlorobiphenyl, 2,2',3,5'-Tetrachlorobiphenyl, 2,2',5,5'-Tetrachlorobiphenyl, 2,3',4',5-Tetrachlorobiphenyl, 2,3',4',6-Pentachlorobiphenyl, 2,3',4,4',5-Pentachlorobiphenyl, 2,2',3,4,4',5'-Hexachlorobiphenyl, 2,2',3,4',5'-Hexachlorobiphenyl, 2,2',4,4',5,5'-Hexachlorobiphenyl, 2,2',4,4',5,5'-Hexachlorobiphenyl, 2,2',3,4,4',5,5'-Hexachlorobiphenyl	Silia <i>Prep</i> C18 Plus 1.5 g / 6 mL SPE-R00830B-06T	Add 0.1 g of L-Ascorbic Acid, 0.35 g of Ethylenediamine Tetra Acetic Acid Trisodium salt and 9.4 g of Potassium Citrate Monobasic to 1 L of surface water.		
Penicillins	Surface water	Amoxicillin, Ampicillin, Penicillin G, Penicillin V, Oxacillin, Cloxacillin, Nafcillin, Dicloxacillin	Silia <i>Prep</i> C18 Plus 200 mg / 3 mL SPE-R00830B-03G	Spike standards in 2 mL of water sample.		
Peptides	Plasma, serum	Angiotensin I, Angiotensin II, Angiotensin III	Silia <i>PrepX</i> HLB 30 mg / 1 mL SPE-P0002-01AA	Add 300 μL of Water to 300 μL of serum. Mix.		
Peptides	Plasma, serum	Oxytocin, Vasopressin	SiliaPrepX WCX 30 mg / 2 mL 96W-P0015-AA	Add 300 μL of 4 % Phosphoric Acid to 300 μL of plasma / serum. Mix.		
Per- and Polyfluoroalkyl Substances (<i>PFASs</i>)	Drinking water	PFBA, PFPeA, PFHxA, PFHpA, PFOA, PFNA, PFDA, PFUdA, PFDoA, PFTrDA, PFTeDA, PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS, FOSA, 4:2 FTS, 6:2 FTS, 8:2 FTS, FOSAA, N-MeFOSAA, N-EtFOSAA, PFHxPA	Silia <i>Prep</i> Amine (<i>WAX</i>) 200 mg / 6 mL SPE-R52030B-06G	Adjust pH between 6.0 and 8.0 with HCl or NaOH (100 mM).		
Phenols	Drinking water	Phenol, 2-Chlorophenol, 2-Methylphenol, 2-Nitrophenol, 2,4-Dimethylphenol, 2,4-Dichlorophenol, 4-Chloro-3-Methylphenol, 2,4-Dirichlorophenol, 2,4-Dinitrophenol, 4-Nitrophenol, 2-Methyl-4,6-Dinitrophenol, 2,4,6-Tribromophenol, Pentachlorophenol	SiliaPrepX DVB 500 mg / 6 mL SPE-P0001-06P	Adjust pH of 1 L water sample to 2.0 with Phosphoric Acid (0.1 M).		
Phthalates	Drinking water	Dimethyl Phthalate, Diethyl Phthalate, Diallyl Phthalate, Dibutyl Phthalate, Diamil Phthalate	Silia <i>Prep</i> C8 500 mg / 3 mL SPE-R31030B-03P	No sample pretreatment required.		

Note: convenient starting points for method development, further optimization may be required



	Chart of Suggested Protocols				
Conditionning	Equilibration	Loading	Washing	Elution	Further Treatment
1 mL of Methanol	1 mL of Methanol / Acetonitrile (50:50)	2 L of water sample Dry the cartridge	No washing required	2 x 0.5 mL of Methanol / Acetonitrile (50:50)	Quantification by GC/MS.
6 mL of Methanol	6 mL of THF 6 mL of n-Pentane	Pretreated sample	3 mL of n-Pentane	2 x 5 mL of Methanol / THF (10:90)	Evaporation under Nitrogen, reconstitution with Methanol and quantification by GC/MS.
3 mL of Methanol	3 mL of Water	Pretreated sample	3 mL of Water Dry the cartridge	1 mL of Methanol	Evaporation under Nitrogen, reconstitution with Water / Methanol (62:38) and identification by HPLC.
5 mL of Ethyl Acetate / Dichloromethane (50:50) Dry the cartridge 10 mL of Methanol	6 mL of Water	Pretreated sample Dry the cartridge	No washing required	5 mL of Ethyl Acetate 5 mL of Dichloromethane	Dry fraction with Sodium Sulfate, evaporation down to 1 mL under Nitrogen and quantification by LC/MS/MS.
4 mL of Methanol	4 mL of Water	Pretreated sample	4 mL of 0.1 % Formic Acid in Water	2 x 2 mL of Acetonitrile	Evaporation under Nitrogen, reconstitution with Water and quantification by LC/MS/MS.
1 mL of Methanol	1 mL of Water	Pretreated sample	1 mL of Methanol / 0.1 % TFA in Water (10:90) Dry the cartridge	1 mL of 2 % Ammonium Hydroxide in Methanol / Water (50:50)	Evaporation under Nitrogen, reconstitution with 0.1 % TFA in Water and identification by HPLC.
1 mL of Methanol	1 mL of Water	Pretreated sample	1 mL of 4 % Phosphoric Acid 1 mL of Methanol / Water (30:70) Dry wells	2 x 0.75 mL of 1 % TFA in Acetonitrile / Methanol (50:50)	Evaporation under Nitrogen, reconstitution with 1 % TFA in 0.1 % Formic Acid / 50 % Methanol in Acetonitrile with 0.1 % Formic Acid (80:20) and quantification by LC/MS/MS.
3 mL of Methanol	5 mL of Phosphate buffer (100 mM, pH 7.0)	Pretreated sample Dry the cartridge	No washing required	6 mL of 1 % Ammonium Hydroxide in Methanol	Evaporation down to 1 mL under Nitrogen and quantification by LC/MS/MS.
6 mL of Dichloromethane	6 mL of Methanol 6 mL of HCI (0.05 N)	1 L of water sample Dry the cartridge	No washing required	6 mL of Dichloromethane	Evaporation down to 1 mL under Nitrogen and identification by HPLC.
3 mL of Methanol	3 mL of Water	200 mL water sample	3 mL of Water	2 x 0.5 mL of Ethyl Acetate	Quantification by GC/MS.



Suggested Protocols for Various Analytes & Matrices (continued)

		Chart of Suggested Proto	cols	
Application	Matrix	Analytes	SPE Cartridge	Sample Pretreatment before Loading
Steroids	Urine, serum	Cortisone, Cortisol, 21-Deoxycortisol, Corticosterone, 11-Deoxycortisol, Fluoxymesterone, Trenbolone, Boldenone, Androstenedione, Nandrolone, Methandienone, 17 α -Hydroxyprogesterone, Testosterone, 16 β -Hydroxystanozolol, Epitestosterone, 5 β -Estran-3 α -ol-17-one, 17 α -Methyltestosterone, Methenolone, 5 α -Estran-3 α -ol-17-one, Norethandrolone, Progesterone, Stanozolol	Silia <i>Prep</i> C8/SAX nec 200 mg / 6 mL SPM-R022830B-06G	Serum: add 4 mL of Phosphate buffer (100 mM, pH 7.0) to 1 mL of sample. Vortex. Urine: add 1 mL of Acetate buffer (100 mM, pH 5.0) and 50 μL of β-glucuronidase to 1 mL of sample. Vortex, heat at 65°C for 1 or 2 hours. Add 2 mL of Phosphate buffer (100 mM, pH 7.0). Vortex.
Sulfonamides	Honey	Sulfanilamide, Sulfathiazole, Sulfamerazine, Sulfamethoxazole, Sulfaquinoxaline	SiliaPrepX SCX 60 mg / 3 mL SPE-P0005-03BB	Add 1 mL of Hydrochloric Acid (2 M) to 1 g of honey, sonicate for 30 minutes and dilute to 5 mL with Citric Acid (300 mM).
Sympathomimetic Amines	Biological fluids (blood, plasma, serum, urine or tissue)	Phenylpropanolamine, Ephedrine, Phentermine, Diethylproprion, Pheniramine, Doxylamine, Chlorpheniramine, Brompheniramine, Diphenhydramine, Amphetamine, Methamphetamine, MDA, MDMA, Pseudoephedrine, Phenylephrine	Silia <i>Prep</i> CleanDRUG 200 mg / 6 mL SPEC-R651230B-06G	Add 3 mL of Phosphate buffer ($100 mM$, $pH 6.0$) to 2 mL of blood / plasma / serum / urine ($or 1 g$ of tissue homogenate). Mix / vortex. Adjust pH to 6.0 ± 0.5 with Sodium Phosphate ($100 mM$). Centrifugate at 2,000 rpm for 10 minutes and discard cellular debris.
Synthetic Dyes & Metabolites	Seafood	Malachite Green, Leucomalachite Green, Crystal Violet, Leucocrystal Violet, Nile Blue, Azure B, Methylene Blue, Brilliant Green, Victoria Blue	Silia <i>Prep</i> CleanDRUG 200 mg / 6 mL SPEC-R651230B-06G	Add 10 mL of 1 % Formic Acid in Acetonitrile and 1 mL of Ascorbic Acid (1 M, used as antioxidant) to 2 g of sample. Vortex for 15 minutes, centrifugate for 10 minutes at 3,000 rpm and 4°C. Take the supernatant and add 20 mL of McIlvaine's buffer (0.1 M, pH 3.5). Vortex and centrifugate again.
Tetracyclines	Honey	Oxytetracycline, Tetracycline, Chlortetracycline	SiliaPrepX WCX 60 mg / 3 mL SPE-P0015-03BB	Dilute honey in Sodium Acetate (100 mM, pH 5.0).
Tobacco	Oral fluid	Cotinine, Anabasine, Nicotine	SiliaPrepX SCX 30 mg / 2 mL 96W-P0005-AA	Collect oral fluid according to instructions on the device, place sample in adequate buffer solution and centrifugate. Add 1 mL of 1 % Formic Acid to 0.5 mL of supernatant.
Topical anesthetics	Serum	Benzocaine, Procaine, Mepivacaine	Silia <i>Prep</i> C8 500 mg / 6 mL SPE-R31030B-06P	Add 500 μL of internal standard solution to 500 μL of serum. Vortex.
Triazoles	Plasma	Fluconazole, Voriconazole, Posaconazole, Ketoconazole, Hydroxy-itraconazole, Itraconazole	Silia <i>PrepX</i> DVB 30 mg / 1 mL SPE-P0001-01AA	Dilute 300 µL of plasma to 1 mL, adjust pH to 2.0.
Urinary Catecholamines	Urine	Dopamine, Norepinephrine, Epinephrine, Normetanephrine, Metanephrine	Silia <i>PrepX</i> WCX 10 mg / 2 mL 96W-P0015-1A	Dilute 75 μL of urine with 150 μL of Ammonium Acetate (<i>250 mM</i>), mix.
Vitamin B3 & Metabolites	Plasma, serum	Nicotinic Acid, Nicotinuric Acid, Niacinamide	SiliaPrep SCX nec 50 mg / 1 mL SPE-R60430B-01B	Add 150 μL of 2 % aqueous Acetic Acid to 50 μL of serum, mix.
Vitamin B7	Serum	Biotin (<i>Vitamin B7</i>)	SiliaPrepX SAX 100 mg / 3 mL SPE-P0010-03C	Add 750 μL of 0.1 % Ammonium Hydroxide in Water to 250 μL of serum. Vortex.

Note: convenient starting points for method development, further optimization may be required



			Chart of Suggested Pro	tocols	
Conditionning	Equilibration	Loading	Washing	Elution	Further Treatment
3 mL of Methanol	3 mL of Phosphate buffer (100 mM, pH 7.0)	Pretreated sample	3 mL of Water 3 mL of Methanol / Water (30:70) Dry the cartridge	2 x 1.5 mL of Methanol	Evaporation under Nitrogen, reconstitution with Methanol / Water (50:50) and quantification by LC/MS/MS.
2 mL of Methanol	2 mL of Water	Pretreated sample	2 x 2 mL of Water 2 x 2 mL of Methanol / Acetonitrile (50:50) Dry the cartridge	2 mL of 2 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water / Acetonitrile (90:10) and quantification by LC/MS/MS.
3 mL of Methanol	3 mL of Water 3 mL of Phosphate buffer (100 mM, pH 6.0)	Pretreated sample	3 mL of Water 3 mL of Acetic Acid (100 mM) 3 mL of Methanol Dry the cartridge	3 mL of Dichloromethane / Isopropanol / Ammonium Hydroxide (78:20:2, pH 11 - 12)	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water and quantification by LC/MS/MS. Or evaporation under Nitrogen, fluoroacylation with PFPA, evaporation under Nitrogen, reconstitution with Ethyl Acetate and quantification by GC/MS.
3 mL of Methanol	3 mL of Water 1 mL of McIlvaine's buffer (0.1 M, pH 3.5)	Pretreated sample (supernatant)	3 mL of 0.1 % Formic Acid in Water 3 mL of 0.1 % Formic Acid in Methanol Dry the cartridge	4 mL of 1 % Triethylamine and 0.5 % Formic Acid in Methanol	Quantification by LC/MS/MS.
2 mL of Methanol	2 mL of Water	4 mL of pretreated sample	4 mL of Water 4 mL of Methanol Dry the cartridge	2 mL of 2 % Formic Acid in Methanol	Evaporation under Nitrogen, reconstitution with Oxalic Acid (10 mM, pH 2.0) / Acetonitrile (80:20) and identification by HPLC.
1 mL of Methanol	1 mL of Water	Pretreated sample	1 mL of Water 1 mL of Acetone / Water (50:50) Dry wells	2 x 500 µL of Ethyl acetate / Isopropanol / Ammonium Hydroxide (70:20:10)	Evaporation under Nitrogen, reconstitution with Ammonium Bicarbonate (20 mM, pH 8.2) / Methanol (90:10) and quantification by LC/MS/MS.
6 mL of Methanol	6 mL of Water	Pretreated sample	6 mL of Water / Methanol (75:25) Dry the cartridge	2 mL of Methanol	Evaporation under Nitrogen, reconstitution with Chloroform and identification by GC.
1 mL of Methanol	1 mL of Water	Pretreated sample	1 mL of 1 % Ammonium Hydroxide in Water 1 mL of Water / Methanol (70:30)	1 mL of Methanol	Evaporation under Nitrogen, reconstitution with Water / Methanol (50:50) and identification by HPLC.
500 μL of Methanol	500 µL of Ammonium Acetate (10 mM)	150 µL of pretreated sample	500 μL of Ammonium Acetate (<i>10 mM</i>) 500 μL of Isopropanol	125 μL of 0.1 % Formic Acid / Isopropanol (85:15)	Quantification by LC/MS/MS.
1 mL of Methanol	1 mL of 2 % Acetic Acid	Pretreated sample	2 x 1 mL of Water / Methanol / Acetic Acid (68:30:2) 2 x 1 mL of 2 % Acetic Acid in Methanol	2 x 400 µL of 5 % Ammonium Hydroxide in Methanol	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid and quantification by LC/MS/MS.
3 mL of Methanol	3 mL of 0.1 % Ammonium Hydroxide in Water	Pretreated sample	3 mL of 0.1 % Ammonium Hydroxide in Water 3 mL of Methanol	2 x 1 mL of 2 % Formic Acid in Methanol	Evaporation under Nitrogen, reconstitution with 0.1 % Formic Acid in Water and quantification by LC/MS/MS.





SPE Accessories

Maximize your Productivity with SiliaPrep Accessories!

SiliCycle offers various accessories for SPE Cartridges and Well Plates to simplify method development and expedite high throughput analysis:

- Vacuum Manifolds
- Empty Tubes
- Adapters & Vacuum Adapters
- 96-Well Collection Plates
- Phase Separator Cartridges
- and Other SPE Accessories...

Silia Prep SPE Vacuum Manifolds

Run multiple samples simultaneously, with a controlled flow rate for higher reproducibility, with Silia*Prep* SPE Vacuum Manifolds. These manifolds are available in 12 and 24-Positions configurations and allow consistent extraction. No possibility of cross-contamination from one sample to another.

The design consists in a clear glass chamber equipped with replaceable individual stopcocks (also known as control valves) and solvent guide needles. The adjustable rack allows the use of a wide variety of collection vessels including 13 and 16 mm test tubes, autosampler vials and volumetric flasks.

Simply apply a vacuum source to elute sample through a cartridge directly to the collection vessel of choice.

Complete sets include:

- Glass chamber, vacuum gauge & bleed valve
- · Cover, gasket, male and female luer fittings
- · Individual stopcocks and needles
- Collection rack with posts, shelves and retaining clips.

SiliaPrep SPE Vacuum Manifolds (Complete Sets)				
Product Number	Description			
AUT-0128-12	12-Positions Silia <i>Prep</i> SPE Vacuum Manifold			
AUT-0129-24	24-Positions Silia <i>Prep</i> SPE Vacuum Manifold			





Silia Prep Vacuum Manifold Accessories

Various replacement parts are available for the two Silia Prep Vacuum Manifolds offered by SiliCycle.

	Silia <i>Prep</i> Vacuum Manifold Accessories					
De	escription		12-Positions Vacuum Manifold		24-Positions Vacuum Manifold	
Sil	ia <i>Prep</i> Vacuum Manifold Complete Set		AUT-0128-12 (1/box)		AUT-0129-24 (1/box)	
Sp	are Parts Ordering Information					
	Glass chamber [Dimensions: Length x Width x Height]		AUT-0182-2 (1/box) [7" x 5.25" x 7"]		AUT-0185 (<i>1/box</i>) [12" x 5.25" x 7"]	
	Vacuum gauge, valve & glass chamber kit	LETE SET	AUT-0187 (1/box)		AUT-0189 (1/box)	
	Top cover, gasket & polypropylene stopcocks kit		AUT-0313 (1/box)		AUT-0315 (1/box)	
Ш	Top cover gasket		AUT-0174 (2/box)		AUT-0193 (2/box)	
E SE	Polypropylene stopcocks		AUT-0146 (12/box)	E SE	AUT-0147 (24/box)	
	Polypropylene needles		AUT-0154 (12/box)	I ⊨	AUT-0155 (24/box)	
COMPLET	Collection rack kit (posts, shelves and retaining clips included)	COMPLET	AUT-0202 (1/box)	COMPLE	AUT-0204 (1/box)	
N C	Plate for 13 mm test tubes	N N	AUT-0205 (1/box)		AUT-0207 (1/box)	
	Plate for 16 mm test tubes		AUT-0208 (1/box)		AUT-0210 (1/box)	
INCLUD	Plate for autosampler vials	INCLUDED	AUT-0213 (1/box)	INCLUD	-	
2	Plate for volumetric flasks		AUT-0214 (1/box)	NC N	-	
	Female luer fittings		AUT-0326 (10/box)		AUT-0326 (10/box)	
	Male luer fittings		AUT-0327 (10/box)		AUT-0327 (10/box)	
	Legs for cover (black)		AUT-0329 (4/box)		AUT-0329 (4/box)	
	Vacuum manifold plugs (yellow)		AUT-0333 (50/box)		AUT-0333 (50/box)	

Note: Stainless Steel needles and Teflon® needles are available upon request.

Silia Prep Waste Containers

Disposable solvent resistant polypropylene containers are available for the 12-Positions manifold. These waste containers greatly simplify sample preparation, solvent disposal and clean-up. Depending on the nature of the solvent used, the waste container can be reused many times prior to discarding.

Note: Waste containers not available for the 24-Positions vacuum manifold.

Silia Prep Drying Manifold Covers

Silia*Prep* Drying Manifold Covers can be used to concentrate samples with a flow of air or gaz (*nitrogen*).

SiliaPrep Drying Manifold Covers			
Product Number	Description		
AUT-0215-12	12-Positions Silia <i>Prep</i> Drying Manifold Cover (<i>1/box</i>)		
AUT-0215-24	24-Positions Silia <i>Prep</i> Drying Manifold Cover (1/box)		



AUT-0176 (10/box)



AUT-0215-12



Silia Prep Adapters

Enable cartridge stacking and easy SPE cartridge connection with syringe or gas lines (for positive pressure).

	SiliaPrep Adapters
Product Number	Description
AUT-0172	SiliaPrep Adapter for 1, 3, 6 & 12 mL SPE (10/box)
AUT-0173	Silia <i>Prep</i> Adapter for 25 & 70 mL SPE (<i>10/box</i>)





AUT-0172 AUT-0173

Silia Prep Vacuum Adapters

Fast, user-friendly, and economical adapters for SPE cartridges. Only a vacuum source is needed.

	SiliaPr	ep Vacuum Adapter - Flasks
Joint	PN	Description
24/40	AUT-0043	24/40 - Silia <i>Prep</i> Vacuum Adapter (<i>1/box</i>)
19/22 AUT-0044		19/22 - Silia <i>Prep</i> Vacuum Adapter (<i>1/box</i>)
14/22 AUT-0045		14/22 - Silia <i>Prep</i> Vacuum Adapter (<i>1/box</i>)

SiliaPrep Vacuum Adapter - Screw Thread Vials						
Thread	PN	Description				
22/400	AUT-0046	22/400 Vial - Silia <i>Prep</i> Vacuum Adapter Without Vial Connector (<i>1/box</i>)				
22/400 AUT-0047		22/400 Vial - Silia <i>Prep</i> Vacuum Adapter With Vial Connector (<i>1/box</i>)				











AUT-0043

AUT-0044

AUT-0045

AUT-0046

AUT-0047

SiliaPrep Empty Tubes

You can use our Silia Prep Empty Tubes to pack your own SPE cartridges with bulk sorbents of your choice.

	Silia <i>Prep</i> Empty Tubes
Formats	Description
SIM-0007-001	Empty 1 mL SPE tube with 2 frits (100/box)
SIM-0008-003	Empty 3 mL SPE tube with 2 frits (100/box)
SIM-0002-006	Empty 6 mL SPE tube with 2 frits (100/box)
SIM-0003-012	Empty 12 mL SPE tube with 2 frits (100/box)
SIM-0004-020	Empty 25 mL SPE tube with 2 frits (100/box)
SIM-0006-060	Empty 60 mL SPE tube with 2 frits (100/box)
SIM-0009-150	Empty 150 mL SPE tube with 2 frits (20/box)




Silia Prep 96-Well Collection Plates


SiliCycle offers Silia*Prep* 96-Well Collection Plates, made from polypropylene with extremely low extractable levels.

These collection plates are available with square deep shape in both 1.0 mL and 2.0 mL well volume (*V-shaped bottom*), and with round bottom in 1 mL only. Cap mats are available for all of these collection plates.


	SiliaPrep 96-Well Collection Plates
Product Number	Description
96W-0009	SiliaPrep 96-Well Collection Plate Square Bottom, 2 mL (50/box)
96W-0010	SiliaPrep 96-Well Collection Plate Square Bottom, 1 mL (50/box)
96W-0011	SiliaPrep 96-Well Collection Plate Round Bottom, 1 mL (50/box)

96-Well Collection Plates Square Shape











2.0 mL Well Volume

96-Well Collection Plates Round Shape





1.0 mL Well Volume

Silia Prep Disposable Reservoir Trays for 96-Well Plates

SiliCycle offers Silia*Prep* Disposable Reservoir Trays to collect waste solvents used during activation, loading and washing steps.


These disposable trays are made of polycarbonate and are compatible with all manifolds used with well plates.


SiliaPrep Disposable Reservoir Trays		
Product Number	Description	
96W-0012	Silia <i>Prep</i> Disposable Reservoir Trays (<i>25/box</i>)	



Silia Prep 96-Well Plate Cap Mats

SiliCycle offers Silia*Prep* 96-Well Plate Cap Mats compatible with most 96-Well Plates available on the market. These cap mats are made of premium-quality silicone, with a PTFE coating for ultra low bleed. Slit and 384-Well Plate cap mats are available under request.







Silia*Prep* 96-Well Plate Square Silicone / PTFE Cap Mats



Silia*Prep* 96-Well Plate Round Silicone / PTFE Cap Mats

	Silia <i>Prep</i> 96-Well Plate Cap Mats				
Well Shape	Quantity Product Number		Description		
	5/box	96M-0001S			
Sauara	25/box	96M-0001S-25	 Silia <i>Prep</i> 96-Well Plate Square Silicone / PTFE Cap Mats		
Square	50/box	96M-0001S-50	(to be used with 96W-0009 & 96W-0010 collection plates)		
	100/box	96M-0001S-100			
	5/box	96M-0001R			
Round	25/box	96M-0001R-25	Silia <i>Prep</i> 96-Well Plate Round Silicone / PTFE Cap Mats		
Roulla	50/box	96M-0001R-50	(to be used with 96W-0011 collection plates)		
	100/box	96M-0001R-100			

Note: Contact us if you are looking for a cap mat not listed above: sampleprep@silicycle.com



Silia Prep Phase Separator Cartridges


SiliCycle offers Silia*Prep* Phase Separator Cartridges to separate the aqueous phase from heavier chlorinated solvents, under gravity. These ready-to-use cartridges are fitted with a proprietary hydrophobic frit and are a great alternative to liquid-liquid extraction, the most popular technique to do this separation. However, this last method is time consuming, requires the use of a glass funnel (*which needs to be washed between each separate extraction*) and is not suitable for multiple extractions. Silia*Prep* Phase Separator Cartridges solve these drawbacks and offer many advantages:

- · Ease of use
- · Efficient and cost saving
- · Compatible with automated systems

Silia <i>Prep</i> Phase Separator Cartridges		
Product Number	Description	
PS-012	SiliaPrep Phase Separator Cartridges, 12 mL (100/box)	
PS-060	SiliaPrep Phase Separator Cartridges, 60 mL (50/box)	
PS-150	Silia <i>Prep</i> Phase Separator Cartridges, 150 mL (<i>25/box</i>)	

Typical Experimental Procedure

- Select the appropriate size of SiliaPrep Phase Separator Cartridge to hold your entire sample volume (both aqueous and chlorinated phases).
- 2. Connect the Silia*Prep* Phase Separator Cartridge on a vacuum manifold. Ensure the collection vessel volume is large enough to entirely recover the organic layer.
 - Note: Do not connect the manifold to a vacuum source
- Transfer the biphasic sample on top of the Silia Prep Phase Separator Cartridge.
- 4. After a few seconds (*under gravity*), the water immiscible chlorinated solvent will start to pass through the frit.
- 5. The proprietary frit used in the Silia*Prep* Phase Separator Cartridge allows the aqueous layer to be left on the column for at least 48 hours without passing through the frit.



Important Advice

· Process under gravity only - Do not apply vacuum or positive pressure

The Silia Prep Phase Separator Cartridges are designed to be used under gravity only. The use of a vacuum or positive pressure source can lead to a loss of separation efficiency.

· Biphasic system required


The sample needs to contain water and a water immiscible solvent (*with greater density than water, to form the lower layer*). Most common solvents are dichloromethane, chloroform and other chlorinated solvents.

Try to minimize the presence of water miscible solvent (*i.e. methanol*, *ethanol* or *acetone*), which can cause problems in obtaining a truly biphasic system. The phase separator may not work effectively if the two phases are merging.

· More efficient compound partition

To obtain a more efficient compound partition between aqueous and organic layers, a liquid-liquid extraction can be done prior to use the phase separator column.





Micro-SPE Tips

Silia Prep Tips for Micro Sample Preparation

- Simple, fast analyte retention & elution with minimal loss
- Sorbents directly embedded into inner cartridge wall
- High binding capacity
- No back-pressure


Overview

Silia Prep Micro-SPE Tips are designed for micro-purification and micro-extraction of femtomole (*fmol*) to picomole (*pmol*) quantities of analytes prior to the analysis by chromatographic techniques and / or mass spectrometry.

The constant improvement in these techniques of analysis has allowed scientists to decrease the limit of quantification in several applications. This lower limit has pushed SPE manufacturers to design new SPE cartridges accepting smaller volumes of analyte.

These tips are specially designed to achieve extraction and purification of small molecules, peptides, phosphopeptides and proteins. They are packed with our Silia*Bond* functionalized silica gels and specialty phases to cover the broadest spectrum of applications requiring small volume of analytes.

The phases are embedded directly in the inner surface of the tip to provide consistent flow rates. Finally, no glue is used during packing procedures in order to prevent any contamination of the analyte.



A - Flow-Through Channel

Eliminates back-pressure.

B - Sorbent

Sorbent is embedded without any glue or binding polymer.

C - Modified Pipette Tip

Easy to fill and release the sample.



SiliaPrep Micro-SPE Tips Sizes

Silia Prep Micro-SPE Tips are available in 3 different cartridge formats, based on the binding capacity of each embedded sorbent.

	SiliaPrep Micro-SPE Tips Specifications				
Tip Volume (μL)	Sample Volume (μL)	Binding Capacity (μg)	Sorbent Weight (μg)	Product Number	
1 - 10	0.5 - 10	1	30	-T1	
10 - 200	2 - 25	2.5	75	-T2	
10 - 200	5 - 50	15	400	-T3	

Silia Prep Micro-SPE Tips are sold in box of 96.



Silia Prep Tips General Experimental Procedure

The following lines present the general experimental procedure for the purification and enrichment of small molecules, peptides and proteins using Silia*Prep* Micro-SPE Tips.

1. Conditionning Step:

Attach the Silia*Prep* Tips to a micropipette. Aspirate / expel the elution solution 5 times and the binding solution 3 times.

2. Loading Step:

Aspirate / expel the sample 20 to 50 times to allow compounds to adsorb onto the sorbent.


3. Washing Step:

Aspirate / expel the binding solution 10 times and discard the expelled solution each time.

4. Elution Step:

Aspirate / expel the elution solution 10 times and collect the expelled solution in a suitable clean tube. Repeat with a fresh portion of elution solution if you want to be sure to collect all of the adsorbed compounds.

Note: repeat 3 - 5 times for carbon black sorbent



Silia Prep Micro-SPE Tips Application

	Micro-Extraction of De	extromethorphan from Plasma		
CARTRIDGE	Silia <i>PrepX</i> Tips C18 10 μL / 30 μg Part Number: SPET-C18-T1			
SAMPLE PRETREATMENT	8 μL of plasma sample was mixed with 2 μL of internal standard (Dextromethorphan-d3 at 10 ng/mL in Methanol)			
CONDITIONNING STEP	8 μL of Methanol (10 aspira	te / expel)		
EQUILIBRATION STEP	8 μL of water (10 aspirate / expel)			
LOADING STEP	Plasma sample (30 aspirate / expel)			
WASHING STEP	8 μL of water (10 aspirate / expel) then 8 μL of 25 % Methanol in water (10 aspirate / expel)			
ELUTION STEP	8 μL of Acetonitrile (30 aspirate / expel)			
FURTHER TREATMENT	Quantification by LDTD/MS/MS (collaboration with Phytronix)			
RECOVERY	at 10 ng/mL			
	Dextromethorphan	86 %		
	Dextromethorphan d-3	80 %		







Silia Prep Tips Sorbent Selection Guide

SiliaPrep Tips Sorbent Selection Guide						
Molecule	Application	Analyte	Sorbent			
	Desalting	All	C18; Carbon Black			
	Protein removal	All	C18; HILIC			
	Metal scavenging	All	Cysteine; DMT; Imidazole; PSA; TAAcOH; TAAcONa; Thiol; Thiourea; Triamine			
		Hydrophobic	C18; HLB; DVB; Carbon Black; HILIC			
Small Molecules		Hydrophilic	Silica; Cyano; Carbon Black; HILIC			
	Enrichment	Neutral	C18; HLB; DVB; Carbon Black; HILIC; Cyano			
	Enrichment	Cationic	SCX; WCX; Polymeric SCX & WCX			
		Anionic	SAX; NH ₂ ; Polymeric SAX & WAX			
		Fluorinated Compounds	Fluoro			
	Desalting	All	C4; C8; C18; Carbon Black; HILIC			
	SDS removal	All	SDS Removal			
Peptides	Enrichment	Glycopeptide	Carbon Black; HILIC; TiO ₂			
		Phosphopeptide	TiO ₂ ; ZrO ₂ ; TiO ₂ /ZrO ₂ ; SAX; NH ₂ ; Polymeric SAX & WAX			
		Other peptide	SAX; NH ₂ ; SCX; WCX; Polymeric SAX, WAX, SCX & WCX			
	SDS removal	All	SDS Removal			
Drotoino	Tryptic digestion	All	Trypsin			
Proteins	Enrichment	Phosphoprotein	TiO ₂ ; ZrO ₂ ; TiO ₂ / ZrO ₂ ; SAX; NH ₂ ; Polymeric SAX & WAX			
	Enrichment	Other protein	SAX; NH ₂ ; SCX; WCX; Polymeric SAX, WAX, SCX & WCX			
	Desalting	All	Carbon Black			
Oligo caccharidas		Sulfated glycan	SAX; XSAX			
Oligo-saccharides	Enrichment	Sialo-glycan	SAX; XSAX			
		Other oligosaccharide	Carbon Black; HILIC; TiO ₂			

Silia Prep XL Tips for Bigger Volumes

For bigger volumes, we also offer Silia*Prep* XL Micro-SPE Tips, in 3 different formats. Please note these tips are top loading instead of by aspiration.

SiliaPrep XL Micro-SPE Tips Specifications				
Tip Volume (μL)	Sample Volume (µL) Binding Capacity (µg) Sort		Sorbent Weight (mg)	Product Number
1 - 10	1 - 10	400	4	-T1
10 - 200	2 - 25	1,000	10	-T2
100 - 1,000	20 - 1,000	5,000	50	-Т3

Silia Prep XL Micro-SPE Tips T1 and T2 are sold in box of 96, T3 in box of 20.



Silia Prep Tips Sorbent Descriptions & Ordering Information

Y).	Silia <i>Prep_</i> Tips S	Sorbent Descriptions and	Ordering Information	n	
				Product Number	
Sorbent	Description		10 μL / 30 μg	200 μL / 75 μg	200 μL / 400 μg
C18	Highest hydrophobic character sorbent. Main and small molecules purification & enrichmen desalting analysis.	SPET-C18-T1	SPET-C18-T2	SPET-C18-T3	
C8	Mid-level hydrophobic sorbent. Mainly used proteins and peptides requiring a lower hydrorotein / peptide desalting analysis.		SPET-C8-T1	SPET-C8-T2	SPET-C8-T3
C4	Lowest hydrophobic character sorbent. Main! & enrichment or protein / peptide desalting ar		SPET-C4-T1	SPET-C4-T2	SPET-C4-T3
HLB	Polymeric sorbent with an hydrophilic-lipophil hydrophobic and neutral molecules enrichme	-	SPET-HLB-T1	SPET-HLB-T2	SPET-HLB-T3
DVB	Higly hydrophobic polymeric sorbent. Mainly neutral molecules enrichment.	used for hydrophobic and	SPET-DVB-T1	SPET-DVB-T2	SPET-DVB-T3
Carbon Black	Hydrophilic and hydrophobic character. Mainl oligosaccharides and other macromolecules of protein / peptide desalting.		SPET-CB-T1	SPET-CB-T2	SPET-CB-T3
HILIC	Moderatly polar sorbent. Mainly used for prote desalting, small molecules enrichment and de spectrum of detergents).		SPET-HI-T1	SPET-HI-T2	SPET-HI-T3
Cyano (CN)	Both polar and hydrophobic character. Mainly neutral molecules enrichment.	used for hydrophilic and	SPET-CN-T1	SPET-CN-T2	SPET-CN-T3
Silica	Most polar sorbent. Mainly used for hydrophil	ic molecules enrichment.	SPET-SI-T1	SPET-SI-T2	SPET-SI-T3
SAX	Strong anion exchanger sorbent. Mainly used	SPET-SAX-T1	SPET-SAX-T2	SPET-SAX-T3	
NH ₂ (WAX)	Weak anion exchanger sorbent. Mainly used (phosphopeptides and phosphoproteins).	SPET-NH2-T1	SPET-NH2-T2	SPET-NH2-T3	
SCX	Strong cation exchanger sorbent. Mainly used	SPET-SCX-T1	SPET-SCX-T2	SPET-SCX-T3	
WCX	Weak cation exchanger sorbent. Mainly used for strong bases enrichment.		SPET-WCX-T1	SPET-WCX-T2	SPET-WCX-T3
SAX Polymeric	Polymeric sorbent functionalized by a strong anion exchanger. Mainly used for weak acids enrichment.		SPET-XSAX-T1	SPET-XSAX-T2	SPET-XSAX-T3
WAX Polymeric	Polymeric sorbent functionalized by a weak anion exchanger. Mainly used for strong acids enrichment (phosphopeptides and phosphoproteins).		SPET-XWAX-T1	SPET-XWAX-T2	SPET-XWAX-T3
SCX Polymeric	Polymeric sorbent functionalized by a strong cation exchanger. Mainly used for weak bases enrichment.		SPET-XSCX-T1	SPET-XSCX-T2	SPET-XSCX-T3
WCX Polymeric	Polymeric sorbent functionalized by a weak cation exchanger. Mainly used for strong bases enrichment.		SPET-XWCX-T1	SPET-XWCX-T2	SPET-XWCX-T3
TiO ₂	High selectivity for multiple phosphorylated per phosphopeptide enrichment and phospholipic	I removal.	SPET-TI-T1	SPET-TI-T2	SPET-TI-T3
ZrO ₂	High selectivity for mono-phosphorylated pep phosphopeptide enrichment and phospholipic		SPET-ZR-T1	SPET-ZR-T2	SPET-ZR-T3
TiO ₂ / ZrO ₂	Excellent alternative for the enrichment of a broad spectrum of phosphopeptides (<i>litterature suggests only 30 % overlap in phosphopeptides isolated by TiO</i> ₂ <i>versus ZrO</i> ₂) and phospholipid removal.		SPET-TIZR-T1	SPET-TIZR-T2	SPET-TIZR-T3
SDS Removal	Used to remove SDS from peptides and proteins.		SPET-SDS-T1	SPET-SDS-T2	SPET-SDS-T3
Trypsin	Used to cleave proteins and peptides at the C-terminal side, with minimal protease contaminants.		SPET-TRYP-T1	SPET-TRYP-T2	SPET-TRYP-T3
Fluoro	Fluorinated sorbent. Mainly used for fluorine coenrichment.	ontaining molecules	SPET-FL-T1	SPET-FL-T2	SPET-FL-T3
		Cysteine	SPET-CYS-T1	SPET-CYS-T2	SPET-CYS-T3
		DMT	SPET-DMT-T1	SPET-DMT-T2	SPET-DMT-T3
		Imidazole	SPET-IMIDAZ-T1	SPET-IMIDAZ-T2	SPET-IMIDAZ-T3
Metal	Choice of 9 metal scavenging sorbents.	Diamine	SPET-PSA-T1	SPET-PSA-T2	SPET-PSA-T3
Scavengers	Mainly used to lower the residual metal concentration of various metal complexes	TAACOH	SPET-TAACOH-T1	SPET-TAACOH-T2	SPET-TAACOH-T3
	(Pd, Pt, Rh, Ru, Ni, Sn, etc).	TAAcONa	SPET-TAACONA-T1	SPET-TAACONA-T2	SPET-TAACONA-T3
		Thiol Thiourea	SPET-THIOLIDEA T1	SPET-THIOLIDEA TO	SPET-THIOLIDEA TO
			SPET-THIOUREA-T1		SPET-THIOUREA-T3
		Triamine	SPET-TRINH2-T1	SPET-TRINH2-T2	SPET-TRINH2-T3

 $\textbf{Note} : \mathsf{Add} \ "XL" \ after \ "SPET" \ for \ ordering \ Silia \\ \textit{Prep} \ XL \ Tips. \ For \ example : SPETXL-C18-T1.$





Sample Preparation & Pesticide Analysis

If you want to improve your sample prep & cleanup procedures, Silia*Fast*™ cartridges offer a simple, economical and highly performant alternative.

FaPEx®: One of the Fastest Extraction / Clean-up Approaches for Pesticide Residue Analysis!



FaPEx stands for "Fast Pesticide Extraction" and may be considered as "QuEChERS made even easier". This 1-step extraction method preceding LC/MS/MS or GC/MS/MS analysis will ensure you:

- Extraction of thousands of pesticides simultaneously
- Reduction of at least 60 % labor cost
- · Faster than existing methods
- Less operating equipment, less organic solvents and waste than QuEChERS
- · Impressive versatility
- · High reliability

Pesticide residues in agricultural products are receiving ongoing and considerable attention from consumers. As global requirements for food safety are growing in numbers and complexity, safer, faster and more refined methods are needed to increase efficiency and overall capacity of pesticide residue testing laboratories. Sample preparation has most frequently been the rate-limiting step of analytic labs, and it is within this context that the FaPEx technology plays its full role.

Currently, the fastest and most internationally recognized rapid extraction technique for pesticide residue analysis is the QuEChERS methodology (see previous page). It is an excellent alternative to more conventional and tedious methods that combine a first Liquid / Liquid Extraction (or LLE) followed by a second Solid Phase Extraction (or SPE).


With overall enhancement in instrument performances (sensitivity, detection ability, resolving power, throughput, other performance attributes), sample preparation techniques had no choice but to keep up with these developments and follow the trend.

A New, Easy, Reliable & Costs Saving Method

FaPEx is an innovative, simple and fast method for the preparation of samples, a crucial step preceding LC/MS/MS or GC/MS/MS analysis.

The FaPEx methodology is based on the same principles than the QuEChERS method, while simplifying it to the extreme. It relies on single-use pre-filled sealed cartridges. After homogenization, the sample is directly treated with a FaPEx cartridge. The resulting cleaned filtrate is ready for injection and chromatographic analysis. Because of its ease of use and great simplicity, FaPEx minimizes handling errors while providing high recovery yields.

More than simply increasing yields of extraction, reducing time and costs associated with pesticide separation, FaPEx also scales down labor time, glassware requirements, specialized equipments and laboratory space requirement. It also significantly diminishes the amount of waste solvents and other chemicals.



Time savings using FaPEx vs standard methods




Design

The basic principle of FaPEx is based on the partitioning of a sample between:

- a liquid phase: the food matrix (*fruits*, *vegetables*, *cereals*, *tobacco*, *etc.*) and key analytes (*pesticides*).
- a stationary phase: the adsorbent (composed of various reagents, depending on the need).

Reagents for removal of matrix interferences:

- dewatering agents: dry the liquid phase and allow a better partition coefficient.
- PSA (Primary Secondary Amine): removes sugars, fatty acids, organic acids, lipids, certain pigments. In combination with a C18 phase, also removes sterols and other types of lipids.
- C18: removes nonpolar organic compounds, saturated hydrocarbons and sterols.
- Carbon black: removes pigments, polyphenols, and other polar organic compounds.







Portfolio

Silia*Fast*[™] is SiliCycle's brand for FaPEx[®] cartridges.

Four different FaPEx products are available according to the various matrices and compound classes below:

	SiliaFast FaPEx Portfolio				
Cartridge	Name	Matrices			
Indicate Process Proce	Silia <i>Fast</i> ™ FaPEx-GEN	General matrices and all forms. Non-pigmented and non-fat containing matrices.			
Southers Street Benefit and Street St	Silia <i>Fast</i> ™ FaPEx-CER	Cereal, rice and grains.			
Southern Ten Proping State of the State of t	Silia <i>Fast</i> ™ FaPEx-PM	Pigmented matrices. Chlorophyll containing matrices.			
STACKES STATE OF STAT	Silia <i>Fast</i> ™ FaPEx-HPM	Highly pigmented matrices. Green tea and other highly pigmented matrices.			

Learn more & Order

Visit www.silicycle.com/fapex

Contact us: sampleprep@silicycle.com



Contact us

SiliCycle Inc.
2500, Parc-Technologique Blvd,
Quebec City (QC) G1P 4S6
CANADA

E-mail: info@silicycle.com

Phone: **+1 418.874.0054**Toll Free **+1 877.745.4292** (North America only)

Fax: **+1 418.874.0355**

Follow us

in www.linkedin.com/company/silicycle

@SiliCycle_Inc

www.facebook.com/silicycle

You Tube www.youtube.com/user/SiliCycle



SiliCycle, a Worldwide Network at your Service



Technical Support

At SiliCycle, we are committed to providing the best technical support possible.

Our worldwide Technical Support Group is comprised of a team of highly qualified M. Sc., Ph. D. Chemists and Engineers, Technical Support Professionals and Service Coordinators who are prepared to troubleshoot, answer questions and provide solutions for your service and applications needs.

In order to better respond to your technical inquiries, feel free to contact us in three different ways:

E-mail: support@silicycle.com Phone: • International +1 418.874.0054

• USA and Canada +1 877.745.4292 (Toll Free)





High Value Silica-based and Specialty Products for Chromatography, Analytical and Organic Chemistry, and Purification

DISCOVER OUR BROCHURES

METAL & ORGANIC SCAVENGING

SiliaMetS® - Metal Scavengers

SiliaBond® - Organic Scavengers

E-PAK® - Fixed Bed Flow-Through Purification Cartridges









CHROMATOGRAPHY & PURIFICATION

SiliaFlash® - Irregular Silica Gels

SiliaSphere™ PC - Spherical Silica Gels

SiliaBond® – Chromatographic Phases

SiliaSep™ - Flash Cartridges

SiliaPlate™ - TLC Plates









SAMPLE PREPARATION

SiliaPrep™ - Silica-based SPE Cartridges & Well Plates

SiliaPrepX[™] - Polymeric SPE Cartridges & Well Plates

SiliaPrep™ - Micro-SPE Tips

SiliaFast™ - FaPEx® Cartridges for Pesticide Residue Analysis









ANALYTICAL & PREPARATIVE CHROMATOGRAPHY

SiliaSphere™ - Spherical Silica Gels

SiliaChrom® - HPLC Columns









ORGANIC SYNTHESIS

SiliaBond® - Reagents & Oxidants









R&D SERVICES

Metal & Organic Scavenging Screening Organic Synthesis Extraction & Purification Custom Column Packing

Material Science













Phone: +1 418.874.0054

Toll Free +1 877.745.4292 (North America only)

