



Spe-ed SCF Helix Line of Supercritical CO₂ and Subcritical Water Systems



**Working with
Supercritical Fluids
the Endless Applications**



**Applied
Separations**

Investigate the Myriad of Applications with the Helix using Existing CO₂ for Supercritical (SC)CO₂ Work

Just as the DNA “double helix” can produce a large variety of different protein configurations, so too does the Applied Separations’ Helix give you the opportunity to make system configurations to investigate the widest variety of supercritical fluid applications.

EXTRACTIONS

- Natural Products
- Foods, Spices, Herbs, Coffee
- Pharmaceuticals
- Fragrances / Essential Oils
- Cannabis / Hemp
- Inorganics: rare earth elements

MAKING NANO PARTICLES

Aerogel drying, critical point drying, hydrogel processing and drying

REACTIONS

- Tissue scaffolding
- Enzymatic Reactions
- Hydrogenation reactions

CLEANING

- Electronic components
- Medical devices
- Machine parts
- 3-D printed parts

IMPREGNATIONS / INFUSIONS

- Textile dyeing
- Medicinals infusions
- Polymers
- Foaming

Liquids extractions/reactions, counter-current, co-current

- Dyeing
- Metal injection molding

PROCESSING

Tissue processing

And now, even more applications using subcritical water

Even doing Accelerated solvent extraction or using gases and liquids other than CO₂


Helix

For Supercritical Fluids AND/OR Subcritical Water



The Helix is made up of several “base” components. The basic components are put together in a variety of standard or custom configurations to make the system perform the applications you need.

With the base system you will be able to use some of the same components to do separations and extractions as well as many other applications. One day do extractions and on another day do something completely different, like infusions or making nano particles.



Basic Requirements for All Extractions

All extractions require just these 3 Applied Separations, Inc. components:

- Base unit (Helix Basic)
- CO₂ Pump
- Pressure vessel assembly

and these utilities:

- Recirculating bath (chiller)
- A source of air delivered at 7 BAR
- Electrical power: 240v or 120v
- Source of liquid CO₂ or water.

The compact Helix Basic is the simplest configuration.

It consists of the:

- CO₂ pump module,
- computer/laptop interface and
- platform holding the pressure vessel

The system can accommodate pressure vessels ranging in size from 24 ml to 1 liter*.

Input, output and vent lines are controlled by shutoff valves located on the front of the unit.

Pressure vessels are heated by specially designed heater assemblies that accurately control the temperature.

An additional CO₂ preheater is employed to ensure that the CO₂ is at the designed temperature before entering the pressure vessel.

A back pressure regulator controls the flowrate of gaseous CO₂

The Helix Basic is controlled through ASI vision software.



The Helix Basic and ALL Helix systems have these features:

- temperatures up to 240°C – no technical limitations because of not having sufficient heat
- pressure up to 10,000 psi (690 BAR) – get maximum density and solvating power, don't leave extracts in the pressure vessel
- holds 24mL to 1L vessels standard; up to 5 liter available – start with a small sample and scale up—no problem
- pump flow rates up to 400mL/min (more with 5 liter vessel add-in) extract quickly and efficiently**
- fully-adjustable, non-clogging, micrometering valves control CO₂ flow rate
- collection into standard glassware, chilled collection to capture volatile compounds, on SPE for further chromatographic separations
- in-vessel chromatographic capabilities
- modifier addition capability – extract polar compounds, too
- multiple flow path – efficiency, efficiency, efficiency
- extract liquid samples as liquids
- stirrer capability –
- view cell...capability to look inside of vessel – keep a video record
- configurations for multiple/varied applications, check applications sections



* larger vessels can also be added in customized systems

**flow rate based on incompressible liquid

www.appliedseparations.com

**Applied
Separations**

Pressure Vessels and Separators

Pressure vessels are heated by specially designed heater assemblies (“clamshell”) that accurately control the temperature.

An additional CO₂ preheater is employed to ensure that the CO₂ is at the designed temperature before entering the pressure vessel.

A back pressure regulator controls the flow rate of gaseous CO₂ if exiting to ambient collection or regulating the pressure in a downstream pressure vessel (e.g. cyclone separator).



Fractionation and Separators

One of the major advantages of supercritical CO₂ is to be able to fractionate samples by manipulating the pressure in the separator. You can go from extracting at a high pressure to discharging an extract at ambient pressure and temperature. Or, the Helix gives you the option of reducing the pressure in a series of steps.

You pick the number of pressure drops - separators - you need for your evaluations. The Helix offers you up to 3 intermediate separators before ending at an ambient condition. With the new Helix configuration, the extractor and the separators are loaded on to a single, operator friendly frame.

**Manipulate the
Supercritical Pressure**

To Change the Density


Which will Change the Solvating Power

**Allowing You to Selectively
Dissolve Compounds**

and Fractionate


The Many Configurations of the Helix

0
CO₂
Separators




Basic Helix
Helix CO₂ Extraction System with
No CO₂ Separators

1
CO₂
Separator




**Helix CO₂ Extraction System with
1 CO₂ Separator**

2
CO₂
Separators



**Helix CO₂ Extraction System with
2 CO₂ Separators**

3
CO₂
Separators



**Helix Extraction System with
3 CO₂ Separators**

Exact configurations of extractors and separators may vary with a system designed for your specific application.


www.appliedseparations.com

**Applied
Separations**

Now let's give you even more options with the Helix... *Subcritical Water*

You can combine both technologies on a single Helix platform. Do both supercritical CO₂ and Subcritical water


0
CO₂
Separators



Helix CO₂ /
Subcritical H₂O Extraction System
with

- No CO₂ Separators
- Dual Use Extractor Vessel
- Subcritical Water


0
CO₂
Separators



Helix CO₂ /
Subcritical H₂O Extraction System
with

- No CO₂ Separators
- Dual Use Extractor Vessel
- Subcritical Water


2
CO₂
Separators



Helix CO₂ /
Subcritical H₂O Extraction System
with

- 2 CO₂ Separators
- Dual Use Extractor Vessel
- Subcritical Water


0
CO₂
Separators



Helix CO₂ /
Subcritical H₂O Extraction System
with Separate CO₂ /H₂O Extractor
Vessels and

- CO₂ Extractor Vessel
- No CO₂ Separator
- Subcritical Water

1
CO₂
Separator



Helix CO₂ /
Subcritical H₂O Extraction System
with Separate CO₂ /H₂O Extractor
Vessels and

- 1 CO₂ Separator
- Subcritical Water

Exact configurations of extractors and separators may vary with a system designed for your specific application.

www.appliedseparations.com

**Applied
Separations**

Helix Basic CO₂ Extraction System with No CO₂ Separator



The Helix Basic give you the opportunity to extract samples up to one liter in volume at up to 10,000 psi and up to 240°C and discharge the extract to an ambient condition.

- temperatures to 240°C – no technical limitations because of not having sufficient heat
- pressure up to 10,000 psi (690 BAR) – get maximum density and solvating power, don't leave extracts in the pressure vessel
- holds 24mL to 1L vessels standard – start with a small sample and scale up—no problem
- pump flow rates up to 400mL/min - extract quickly and efficiently*
- fully-adjustable, non-clogging, micrometering valves control flow rate
- collection into standard glassware, chilled collection to capture volatile compounds, on SPE for further chromatographic separations
- in-vessel chromatographic capabilities
- modifier addition capability – extract polar compounds, too
- multiple flow path – efficiency, efficiency, efficiency
- extract liquid samples as liquids
- stirrer capability
- view cell...capability to look inside of vessel – keep a video record
- configurations for multiple/varied applications, check applications sections



Extractor Module

Pump Module



Ambient
CO₂
Collection



One of the many thousands of applications might be the collection of soluble oils from natural products. Soluble oils may be dissolved in high pressure supercritical CO₂ and collected when depressurized to atmospheric conditions. For example, soluble lipid components from rosehip seeds, flax seeds etc can be easily dissolved in CO₂ and collected in an atmospheric collection container.

*flow rate based on incompressible liquid

www.appliedseparations.com


**Applied
Separations**

Helix CO₂ Extraction System with 1 CO₂ Separator

The Helix CO₂ Extraction System with 1 CO₂ separator give you the opportunity to extract samples up to one liter in volume at up to 10,000 psi and up to 240°C and discharge the extract to a separator or an ambient condition.



- temperatures to 240°C – no technical limitations because of not having sufficient heat
- pressure up to 10,000 psi (690 BAR) – get maximum density and solvating power, don't leave extracts in the pressure vessel
- holds 24mL to 1L vessels standard; – start with a small sample and scale up— no problem
- pump flow rates up to 400mL/min - extract quickly and efficiently*
- fully-adjustable, non-clogging, micrometering valves control flow rate
- collection into standard glassware, chilled collection to capture volatile compounds, on SPE for further chromatographic separations
- in-vessel chromatographic capabilities
- modifier addition capability – extract polar compounds, too
- multiple flow path – efficiency, efficiency, efficiency
- extract liquid samples as liquids
- stirrer capability
- view cell - capability to look inside of vessel – keep a video record
- configurations for multiple/varied applications, check applications sections

An example might be cloves contain 2 major components that may be dissolved in supercritical CO₂ and separated. A dark oily resin can be collected in separator 1 and an essential oil fraction in the final atmospheric collector.



Separator/Extractor Module

Pump Module



Exact configurations of extractors and separators may vary with a system designed for your specific application.

*flow rate based on incompressible liquid

www.appliedseparations.com



Helix CO₂ Extraction System with 2 CO₂ Separators

The Helix CO₂ Extraction System with 2 CO₂ separators gives you the opportunity to extract samples normally up to one liter in volume at up to 10,000 psi and up to 240°C and discharge the extract to 2 different separators; cascading in decreasing pressure from separator #1 to separator #2 and then to ambient collection.

- temperatures to 240°C – no technical limitations because of not having sufficient heat
- pressure up to 10,000 psi (690 BAR) – get maximum density and solvating power, don't leave extracts in the pressure vessel
- holds 24mL to 1L vessels standard;— start with a small sample and scale up—no problem
- pump flow rates up to 400mL/min - extract quickly and efficiently*
- fully-adjustable, non-clogging, micrometering valves control flow rate
- collection into standard glassware, chilled collection to capture volatile compounds, on SPE for further chromatographic separations
- in-vessel chromatographic capabilities
- modifier addition capability – extract polar compounds, too
- multiple flow path – efficiency, efficiency, efficiency
- extract liquid samples as liquids
- stirrer capability
- view cell - capability to look inside of vessel – keep a video record
- configurations for multiple/varied applications, check applications sections

Where this configuration might be used is the extraction black pepper. Black Pepper contains several distinct components that may be dissolved in supercritical CO₂ at 500 bar and 90°C and separated as the pressure is sequentially lowered.



- A resinous material is collected in separator #1,
- a piperine concentrate in separator #2 and
- an essential oil in the final collection container at atmospheric pressure.



Exact configurations of extractors and separators may vary with a system designed for your specific application.

Separator/Extractor Module

Pump Module



*flow rate based on incompressible liquid

www.appliedseparations.com

Applied Separations


Helix Extraction System with 3 CO₂ Separators

The Helix CO₂ Extraction System with 3 CO₂ separators gives you the opportunity to extract samples up to one liter in volume at up to 10,000 psi and up to 240°C and discharge the extract to 3 different separators; cascading in decreasing pressure from separator #1 to separator #2, to separator #3 and then to ambient collection.

- temperatures to 240°C – no technical limitations because of not having sufficient heat
- pressure up to 10,000 psi (690 BAR) – get maximum density and solvating power, don't leave extracts in the pressure vessel
- holds 24mL to 1L vessels standard;– start with a small sample and scale up—no problem
- pump flow rates up to 400mL/min - extract quickly and efficiently*
- fully-adjustable, non-clogging, micrometering valves control flow rate
- collection into standard glassware, chilled collection to capture volatile compounds, on SPE for further chromatographic separations
- in-vessel chromatographic capabilities
- modifier addition capability – extract polar compounds, too
- multiple flow path – efficiency, efficiency, efficiency
- extract liquid samples as liquids
- stirrer capability
- view cell - capability to look inside of vessel – keep a video record
- configurations for multiple/varied applications, check applications sections

Many CO₂ soluble polymers may be fractionated by molecular weight including silicones, polycarbosilanes, and perfluoropolyethers. For example, polydimethylsiloxane with an average molecular weight of approximately 90,000 can be dissolved in supercritical CO₂ at 80°C and 600 bar and fractionated as the pressure is reduced sequentially from 450 bar to atmosphere.




- Separator# 1 contains an average molecular weight of 150,000,
- Separator #2 contains an average molecular weight of 75,000
- Separator#3 contains an average molecular weight of 46,000
- The ambient collector contains an average molecular weight of 1500.



Separator/Extractor Module

Exact configurations of extractors and separators may vary with a system designed for your specific application.

Pump Module



*flow rate based on incompressible liquid

www.appliedseparations.com

Applied Separations

Helix CO₂ /Subcritical H₂O Dual Use Extraction System with No CO₂ Separators


In this Helix configuration with SCCO₂ and Subcritical H₂O, the same extractor vessel is used for both types of extractions. They cannot be used simultaneously. However, after clean-up from one operation, the other technique can be employed.

The supercritical CO₂ part of the configuration is the same as the Basic Helix, shown on page 8.

In order to do subcritical water extractions, the Helix employs:

- a separate liquid pump to pump the water into the extractor vessel
- a separate collection/condensation system to safely condense the “superheated steam” and collect the desired extracts
- Applied Separations, Inc. standard vessels configured to subcritical water achieve 240°C. Optionally, vessels can be configured to achieve 300°C.

A potential operation would be to perform a supercritical CO₂ extraction like doing a cleanup and then do a subcritical water extraction like extracting flavonoids and polyphenolics from orange peels including the flavonones, hesperidin and narirutin.



Extractor Module



Exact configurations of extractors and separators may vary with a system designed for your specific application.

Helix CO₂ /Subcritical H₂O Dual Use Extraction System with 1 CO₂ Separator

In this Helix configuration with SCCO₂ and Subcritical H₂O, the same extractor vessel is used for both types of extractions. They cannot be used simultaneously. However, after clean-up from one operation, the other technique can be employed.

The supercritical CO₂ part of the configuration is the same as the Helix with one separator, shown on page 9.


In order to do subcritical water extractions, the Helix employs:

- a separate liquid pump to pump the water into the extractor vessel
- a separate collection/condensation system to safely condense the “superheated steam” and collect the desired extracts
- Applied Separations, Inc. standard vessels configured to subcritical water achieve 240°C. Optionally, vessels can be configured to achieve 300°C.

A potential operation would be to perform a supercritical CO₂ extraction like doing a cleanup followed by a simultaneous subcritical water extraction of phenolics (eriocitrin) and essential oils (menthol) from peppermint leaves.



Separator/Extractor Module



Exact configurations of extractors and separators may vary with a system designed for your specific application.

Helix CO₂ /Subcritical H₂O Dual Use Extraction System with 2 CO₂ Separators

In this Helix configuration with SCCO₂ and Subcritical H₂O, the same extractor vessel is used for both types of extractions. They cannot be used simultaneously. However, after clean-up from one operation, the other technique can be employed.


The supercritical CO₂ part of the configuration is the same as the Helix with two separators, shown on page 8.

In order to do subcritical water extractions, the Helix employs:

- a separate liquid pump to pump the water into the extractor vessel
- a separate collection/condensation system to safely condense the “superheated steam” and collect the desired extracts
- Applied Separations, Inc. standard vessels configured to subcritical water achieve 240°C. Optionally, vessels can be configured to achieve 300°C.



Separator/Extractor Module



Exact configurations of extractors and separators may vary with a system designed for your specific application.

Helix CO₂ /Subcritical H₂O Extraction System with

- **Separate Dedicated CO₂ Extractor Vessel**
- **No CO₂ Separator**
- **Separate Dedicated H₂O Extractor Vessels**

In this Helix configuration with SCCO₂ and Subcritical H₂O, there are two extractor vessels. One is dedicated to supercritical CO₂ and the other is used for subcritical water. Although requiring attention, it is possible to operate both techniques simultaneously.

The supercritical CO₂ part of the configuration is the same as the Basic Helix, shown on page 8.

In order to do subcritical water extractions, the Helix employs:



- a separate liquid pump to pump the water into the extractor vessel
- a separate collection/condensation system to safely condense the “superheated steam” and collect the desired extracts
- Applied Separations, Inc. standard vessels configured to subcritical water achieve 240°C. Optionally, vessels can be configured to achieve 300°C.

A potential operation would be to perform a simultaneous supercritical CO₂ extraction while also doing a subcritical extraction of polyphenols and flavonoids from grape pomace.



Separator/Extractor Module

Pump Module



Condenser

Ambient CO₂ Collection



H₂O Pump



H₂O Collection



Exact configurations of extractors and separators may vary with a system designed for your specific application.

Helix CO₂ /Subcritical H₂O Extraction System with

- **Dedicated CO₂ Extractor Vessel**
- **Dedicated H₂O Extractor Vessel**
- **1 CO₂ Separator**

In this Helix configuration with SCCO₂ and Subcritical H₂O, there are three extractor vessels. One is dedicated to supercritical CO₂, one is used for subcritical water, and the other is a CO₂ separator. Although requiring attention, it is possible to operate both techniques simultaneously.

The supercritical CO₂ part of the configuration is the same as the Basic Helix, shown on page 8.

In order to do subcritical water extractions, the Helix employs:

- a separate liquid pump to pump the water into the extractor vessel
- a separate collection/condensation system to safely condense the “superheated steam” and collect the desired extracts
- Applied Separations, Inc. standard vessels configured to subcritical water achieve 240°C. Optionally, vessels can be configured to achieve 300°C.


A potential operation would be to perform a simultaneous supercritical CO₂ extraction while also doing a subcritical extraction of

- theobromine, caffeine, theophylline,
- epicatechin, catechin, chlorogenic acid and gallic acid from waste cocoa shells including
- valuable degradation products such as 5-hydroxymethylfurfural (5-HMF), furfural, levulinic acid, lactic acid and formic acid.



Separator/Extractor Module

Exact configurations of extractors and separators may vary with a system designed for your specific application.




Why Supercritical Fluids?

Carbon dioxide is in its supercritical fluid state when both the temperature and pressure equal or exceed the critical point of 31°C and 73 atm (see diagram). In its supercritical state, CO₂ has both gas-like and liquid-like qualities, and it is this dual characteristic of supercritical fluids that provides the ideal conditions for extracting compounds with a high degree of recovery in a short period of time.

Supercritical Fluids “Green” Revolutionize Your Processes

No longer an exotic laboratory curiosity, but now a cost-effective tool to improve your process development.

Use “green” existing carbon dioxide and no solvents



By controlling or regulating pressure and temperature, the density, or solvent strength, of supercritical fluids can be altered to simulate organic solvents ranging from chloroform to methylene chloride to hexane.

This dissolving power can be applied to purify, extract, fractionate, infuse, and recrystallize a wide array of materials.

Because CO₂ is non-polar, a polar organic co-solvent (or modifier) can be added to the supercritical fluid for processing polar compounds. By controlling the level of pressure/temperature/modifier, supercritical CO₂ can dissolve a broad range of compounds, both polar and non-polar.


Not only can the supercritical fluid act as a solvent, it also has the characteristics of a gas. One can apply changes in pressure to effective foaming, diffusion, and many other additional physical parameters. The applications are endless, limited only to one's imagination.

Why Subcritical Water?

Water can be physically altered to dissolve non polar compounds

There is a great demand for environmentally friendly extraction procedures that eliminate the need for toxic organic solvents. Supercritical carbon dioxide is one and has been highlighted here. However, one can also use water in some innovative ways.

Water can be held in a liquid state above its normal boiling temperature (100°C) by increasing the pressure. This process is demonstrated using a simple pressure cooker. Water above its normal boiling point but below its supercritical temperature (374°C), is “subcritical” water.



Water can be used as an alternative to toxic organic solvents because the dissolving power of water can be changed by changing its temperature. Water can be physically altered to become “subcritical” water

Polar

- Water at room temperature has a high dielectric constant - dissolving polar and ionic compounds.

More Non Polar

- Water at high temperatures has a lower dielectric constant - dissolving non-polar organic compounds

Viscosity and surface tension decrease

- Viscosity and surface tension decrease with increasing temperature
- Solubility of solutes increases
- Solute diffusion increases, decreasing time to solubilize

Replace alcohol with subcritical water

- Since the dielectric constant of water decreases significantly as it is heated, it can behave like alcohol

Subcritical water may be used to extract many organic molecules

- Eliminate the use of conventional organic solvents.

Variation of Dielectric Constant with Temperature for H₂O

