

Counting with EVE, affordable automatic cell counting

Benchtop automatic cell counter, *EVE*! No more manual counting!



EVE, Automatic cell counter, with states-of-the-art optics provides image analysis for cell counting. The EVETM is a benchtop size counter, designed to measure live, dead, and total counting, which gives us capability for viability with accuracy and precision.

O Key features and benefits

Fine distinction of clumped cell

· An advanced analysis algorithm provides accurate results in clumped cell

Applicable for broad range of cell sizes and types

• From primary (tissue, blood) cell lines to stem cells

User friendly

• 7 inch LCD touch screen • Benchtop size • No maintenance needed

Data storage and analysis

* Stores 500 test results, transfers the data to PC using USB drive * Provided PC software will give you more in-depth analysis


Cell size gating function

· Select range of cell size


O Convenient for use



O Correlation of EVE™ and manual counting




Measuring from the EVE™ extends further the high concentration range than hemocytometer readings




More data window to check the detail result and select the cell size gating

Faster! Accurate! Easy to use!

o Total cell counting result (Clumped cell)



The MCF7 breast cancer cell line can be clumpy. The EVE™ counting algorithm identifies and counts individual cells within these cell clumps for accurate analysis.



Clumped cells (HepG2, MCF7, F9) were counted with EVE[™], a hemocytometer and the competitor's automated cell counter. In accuracy and precision, EVE[™] is comparable to the hemocytometer for all cell lines. The competitor's cell counter was shown significantly less numbers in clumped cell. Cell counting results were obtained seven times repetition.

o Cell lines validated on EVE™

Adherent
Adherent
Adherent
Suspension
Suspension
Adherent

Specifications

Item	Description
Countkig time	< 20seconds
Cell measurement range (cells/mL)	1x10 ⁴ -1x10 ⁷
Optimal measurement range (cells/mL)	1x10 ⁵ -4x10 ⁶
Cell size range	5-60 μm
Sample volume	10 μL
Staining Method	Trypan blue stain
Display	7" LCD touch screen
Exporter fomats	JPEG (image), CSV(raw data)
Data export	USB drive
Dimensions	27cm(w)x20cm(d)x19cm(h)
Weight	2.1 kg (4.6 lb)

Ordering information

Cat No.	Description	Contents
EVE-MC	EVE™ Automatic cell counter	Main device, Power supply, USB drive, EVE™ Cell counting slide (50 slides)
EVS-050	EVE™ Cell counting slide	50 slides (100 counts)
EVS-1000	EVE™ Cell counting slide	1,000 slides (2,000 counts)
EVS-5000	EVE™ Cell counting slide	5,000 slides (10,000 counts)
EBB-001	Test beads Concentration (avg.) 1.0 x 10 ⁶	1 mL





Scan for more information.

