



Analysis of N-Nitrosamines and Nitrosatable Substances From Elastomer or Rubber Teats, Soothers, and Toys

AS PER BS EN 71-12:2016 AND BS EN 12868:2017

Testing for Nitrosamines

Nitrosamines are a group of carcinogenic and mutagenic compounds. Daily and prolonged contact with manufactured items such as toys and soothers, will lead to continual exposure of these compounds. Potentially any item containing nitrites/nitrates could contain nitrosamines and any process that uses secondary amines, where the manufacturing process results in an introduction of atmospheric nitrous oxide, could produce nitrosamines.

Concern over cumulative exposure has resulted in Ellutia developing an analysis solution, with specific focus on exposure to babies and young children under the scope of UK method BS EN 71-12:2016 (Safety of toys; N-nitrosamines and N-nitrosatable substances). This method utilises an Ellutia 200 GC for separation followed by an Ellutia 810 TEA (Thermal Energy Analyser) for detection. TEA is a BS EN 12868:2017 approved method for analysing compounds containing nitroso groups.

In accordance to 93/11/ECC, total concentration for N-nitrosamines or N-nitrosatable substances must not exceed 10 µg/kg individually or 100 µg/kg cumulatively, furthermore, current BS EN 71-12:2016 emphasises that baby products should contain less than 10 µg/kg of any individual N-nitrosamine and less than 20 µg/kg of total N-nitrosamines per item sold.

Method

The BS EN 71-12:2016 method covers 13 different nitrosamines including some low or non-volatile nitrosamines such as N-nitrosodiethylamine (NDELA). This can be overcome by volatilisation using a commercially available derivatisation kit containing BSTFA + 1% TMCS.

GC Conditions

Injector Temperature:	220°C
Linear Type:	Focus Liner with Wool
Carrier Gas Type:	Helium
Carrier Gas Control Method:	Constant Flow
Gas Flow:	2psi
Injection Type:	Split-less
Split-less Time:	0.5min
Column Type:	EL-FFAP 15m x 0.32mm x 0.5 µm

Detector Conditions

Detector:	810 TEA
Detector Temperature:	250°C
Pyroliser Temperature:	500°C
Detector Sensitivity:	230
Vacuum:	0.39 torr
Oxygen Flow:	2psi

Column Temperature Program

Initial Temperature:	60°C (hold 0 min)
Temperature Ramp 1:	10°C min ⁻¹ to 160°C (hold 0 min)
Temperature Ramp 2:	20°C min ⁻¹ to 160°C (hold 3 min)

Results

A variety of children and baby items that would fall within the scope of BS EN 71-12 were sampled and extracted using the Ellutia 810 TEA (Thermal Energy Analyser) and quantitation limits were calculated to prove the instrumentation achieves the necessary Limits of Quantitation (LOQ) in accordance with BS EN 12868:2017. The results shown in Table 1 and Figure 1 show how nitrosamines can be successfully separated by the 200 Series GC and detected with the 810 Series TEA with easily matching or exceeding Limits of Quantitation stated in the EN 12868:2017 method.

Part No	Name	Abbreviation	CAS Number	BS EN 1268:2017 LOQ ($\mu\text{g}/\text{kg}$)
1	N-nitrosodimethylamine	NDMA	62-75-9	1.0
2	N-methylethylamine	NMEA	624-78-2	1.0
3	N-nitrosodiethanamine	NDEA	55-18-5	1.0
4	N-nitrosodipropylamine	NDPA	621-64-7	1.0
5	N-nitrosodipropylamine	NDBA	924-16-3	1.0
6	N-nitrosopiperidine	NPIP	100-75-4	1.0
7	N-nitrosopyrrolidine	NPYR	930-55-2	1.0
8	N-nitrosomorpholine	NMOR	59-89-2/67587-56-8	1.0
9	N-nitroso N-ethyl N-phenylamine	NEPhA	612-64-6	5.0
10	N-nitroso N-methyl N-phenylamine	NMPhA	614-00-6	5.0
11	N-nitrosodibenzylamine	NDBzA	5336-53-8	5.0

Table 1: List of N-nitrosamines used with subsequent BS EN 12868:2017 Quantitation Limit requirements.




Figure 1: A chromatogram showing the separation of a multi-component nitrosamine standard.

Concentration measured in $\mu\text{g}/\text{kg}$			
	NDELA	NDMA	NDEA
Blank IS	(2.1)	(3.4)	-
Finger Paints	3.3	2.2	-
Toy 1	0.7	0.3	-
Toy 2	0.8	0.1	-
Looms	1.4	0.2	-
Soothers	0.8	0.5	-
Balloons	-	1.5	0.2

Table 3: Concentration of Nitrosamines found in toys and soothers ($\mu\text{g}/\text{kg}$).

The legislation also requires the use of a deuterated internal standard, however as this analysis is GC-TEA and not GC-MS, the internal standards used in this study were not necessary. A blank reading was taken for the solutions with an internal standard. Spike and recovery calculation were used to determine the analyte detection, as seen in Table 3.

Current EU and U.S. legislations (FDA and CPSC) for nitrosamines recommend that baby products should contain less than $10 \mu\text{g}/\text{kg}$ of any individual and $20 \mu\text{g}/\text{kg}$ of total nitrosamines per item sold. The above data presents that, whilst most products tested showed the presence of nitrosamines, N-nitrosodimethylamine (NDMA) and NDELA were observed in every tested sample at low levels, whereas N-nitrosodiethylamine (NDEA) was only found in the balloons. Therefore, none of the samples tested were found to contain levels above the legislation for individual and total nitrosamines.



Conclusion

As the TEA is a Nitrogen specific detector, the sensitivity and selectivity for the compounds tested is superior to other detectors. The ease of set up and use sets the TEA aside from more costly and complicated systems.


Separation of N-nitrosamines within the BS EN 12868:2017 recommended LOQs can be performed on Ellutia 810 Thermal Energy Analyser. Baseline separation and reproducible chromatography was obtained for the compounds within the concentration range of $1 - 5 \mu\text{g}/\text{kg}$ with LODs and LOQ matching or exceeding those stated in EN 12868:2017.

Further monitoring of nitrosamines within toys and other rubber/latex products is essential. Whilst manufacturers continue to reduce the exposure risk, it is paramount to reduce baby's and children's exposure to these carcinogenic compounds through strict monitoring and control. The Ellutia 810 Thermal Energy Analyser system manifests compliance with appropriate European legislations and the method demonstrates enhanced sensitivity in ability to detect and quantify these harmful compounds.

Ordering Guide

Main Instruments	
TEA 810 - 115V or 230V	(Part no. 32000810)
Cold Trap Replacement	(Part no. 32090001)
200 Series GC - 110V	(Part no. 20501130)
200 Series GC - 230V	(Part no. 20500130)
Elution Software, Single Instrument	(Part no. 23001001)
EL-FFAP 15m x 0.32mm x 0.5µm Capillary Column	(Part no. 51100343)
Optional	
Ellutia EL3100A Automatic Liquid Autosampler - 15 position	(Part no. 30500011)
Ellutia EL3000A Automatic Liquid Sampler - 121 position	(Part no. 30500010)
Ellutia EL3200A Automatic Liquid Sampler - 209 position	(Part no. 30500012)
Autosampler Control Software	(Part no. 23001012)
Accessories	
2ml Vials	(Part no. 20511101)
2ml Vial Screw Caps	(Part no. 20511107)
1 µl Liquid Syringe	(Part no. 20511204)
Septa	(Part no. 51100298)
20ml Headspace Vial and Cap	(Part no. 30501003)
Solutions Required	
BSTFA	
TMCS	





To learn more about Ellutia's range of chromatography solutions, please visit:

www.ellutia.com




Colston House, 200 Lancaster Way Business Park, Ely, Cambridgeshire, CB6 3NX, UK

Tel: +44 (0)1353 669916 Web: www.ellutia.com

Ellutia Limited Registered in England Number 2967460

Registered Address Colston House, 200 Lancaster Way Business Park, Ely, Cambridgeshire, England, CB6 3NX



SUBSCRIBE



CONNECT



SOLVE

