

# mini vortex

## MIXER

### SMALL SIZE, POWERFUL MIXING

With touch-activated operation, the Mini Vortex Mixer provides powerful mixing in a small footprint for tubes up to 50 mL.

### SMALL SIZE, POWERFUL MIXING

The Mini Vortex Mixer delivers a fixed speed of 2,800 rpm by only applying slight pressure to the center of the spill-resistant vortex head. Remove the pressure and the vortex motor immediately stops. The entire operation is quiet and smooth and contained in a small 10 x 10 cm (4 x 4 in.) footprint.



#### MINI VORTEX MIXER

##### Specifications

Specification	Description
Speed	Fixed, 2800 rpm
Operating Mode	TOUCH mix
Head	Rubber cup for single tube mixing
Orbit	4 mm
Operating Environment	4°C to 65°C
Dimensions (WxDxH)	9.4 x 9.9 x 6.6 cm / 3.7 x 3.9 x 2.6 in.
Weight	0.4 kg / 1 lb.
Electrical	universal, 100V-240V, 50-60 Hz

##### Ordering Information

Description	Part Number
Mini Vortex Mixer, 100-240V, uS Plug	36110700
Mini Vortex Mixer, 100-240V, Eu plug	36110710
Mini Vortex Mixer, 100-240V, uK plug	36110720