

roto-mini PLUS

FLEXIBLE VARIABLE SPEED ROTATOR FOR A RANGE OF APPLICATIONS

With speeds from 5 to 70 rpm, the Roto-Mini Plus variable speed rotator offers three operational mixing modes (Mix, Pause, and Rock) and accommodates standard tube sizes from 0.5 mL to 50 mL, making it a flexible addition to any lab.

FLEXIBLE OPERATIONAL MODES

With speeds from 5 to 70 rpm, the Roto-Mini Plus variable speed rotator provides three different mixing modes to cover many applications, including:

- Mix – samples rotate with oscillations for increased mixing
- Pause – sample rotation includes periodic pauses
- Rock – samples gently rock in a see-saw manner.

ADAPT TO MANY TUBE STYLES

The Roto-Mini Plus includes four holders that cover the following common tube sizes:

- 11 x 1.5/2.0 mL tubes
- 5 x 5/7 mL tubes
- 5 x 10/15 mL tubes
- 6 x 50 mL tubes/24 x 0.5 mL tubes



ROTO-MINI PLUS

Description	Part Number
Roto-Mini Plus with tube holders, 115V, u S Plug	36110100
Roto-Mini Plus with tube holders, 230V, Eu Plug	36110110
Roto-Mini Plus with tube holders, 230V, u K plug	36110120