



vortex

MIXER

MAXIMIZE MIXING, MINIMIZE NOISE AND VIBRATION

With a unique balance system, our powerful Vortex Mixer offers instant vortex mixing without the noise and vibration. The mixer features both touch and continuous operation, variable speeds, and can be used in a range of temperature environments.

MAXIMIZE MIXING, MINIMIZE NOISE

The Vortex Mixer relies on a novel counter-balance system that provides powerful mixing action with decreased noise and vibration over other mixing instruments, even for 50 mL tubes. With its suction feet, the Vortex Mixer remains stable on your bench.

The Vortex Mixer offers instant variable speed control from 200 to 3,200 rpm and can be used in touch or continuous operating modes.

VERSATILE

The Vortex Mixer operates in a wide range of temperatures (+4°C to 65°C) and features spill-proof electronic controls. Several optional heads featuring both flat and multi-tube configurations are also available.



ACCESSORIES

Description	Part Number
Vortexer Flat head	36117310
COMBO™ Head	36117320
Optional Foam Rack for 12 x 15mL tubes	36117330
Optional Foam Rack for 5 x 50 mL tubes	36117340
Horizontal Head, for 12 x 1.5 mL	36117350
Horizontal Head, for 4 x 15 mL	36117360
Horizontal Head, for 2 x 50 mL	36117370

VORTEX MIXER

Description	Part Number
Vortex Mixer with standard cup head, 115V, uS Plug	36110740
Vortex Mixer with standard cup head, 230V, Eu Plug	36110750
Vortex Mixer with standard cup head, 230V, uK Plug	36110760