

# MyPipetman®

SELECT

## IMPROVED PIPETTING AND PERSONALIZATION OPTION

Air-Displacement Pipettes



**Ergonomic, safe, and reliable, MyPIPETMAN® Select offers volume setting choice and optional handle personalization for optimal use in various protocols.**

### Select, Set, and Secure Volumes with the Trilock™ System

The unique and patented Trilock™ volume-locking system enables you to choose different positions when setting the volume:

- Free volume setting position for quick and effortless volume adjustment
- Light ratchet volume setting for precise volume adjustment
- Volume locked position for secure pipetting

With this new Trilock system, the volume is adjusted simply by turning the push button of the pipette.



### Ergonomic Design

Lightweight and well-balanced, MyPIPETMAN Select has an elegant design with a textured effect for a better grip.

With no compromise on quality and performance, the built-in ejection system provides reduced tip ejection forces. This patented ejection system reinforces the pipette's ergonomics for convenient use throughout the day.

### Handle Printing with Your Text on Demand

MyPIPETMAN Select can have their handles printed on demand. Identify your pipette to avoid any risk of contamination by dedicating it to a specific liquid or workflow step and improve your lab Lean Management process.





## Extra Personalization and Accessories

Other accessories like colored windows, tip holders, and tip ejectors are available in different colors and material for additional customization and ease of use.

Standard accessories, such as pipette stands and holders, are also available. Please visit [www.gilson.com](http://www.gilson.com) for ordering information.



Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.



PIPETMAN DIAMOND® Tips are designed to fit your MyPIPETMAN perfectly and guarantee maximum accuracy and precision.

Benefit from a special care with our flexible, tiered service level offerings, tailored to fit your lab's needs, while increasing the lifespan of your pipettes.

MyPIPETMAN is also available for Enterprises offering, both functional and product customization options. Benefit from our Gilson brand and quality standards to build your own specific line of pipettes.

Please visit [www.gilson.com](http://www.gilson.com) for more information on MyPIPETMAN Enterprise

Colored Windows	
Model	Part Number
Window 3digit Dark Blue MyPIPETMAN	FP070101
Window 3digit White MyPIPETMAN	FP070102
Window 3digit Blue MyPIPETMAN	FP070103
Window 3digit Red MyPIPETMAN	FP070104
Window 3digit Yellow MyPIPETMAN	FP070105
Window 3digit Green MyPIPETMAN	FP070106

Accessories	
Description	Part Number
CARROUSEL™ pipette Stand (7 pipettes)	F161401
TRIO™ Pipette Stand (3 pipettes)	F161405
Universal Multichannel Stand	F161417
SINGLE™ Pipette Holder	F161406

MyPIPETMAN® Select				
Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (µL)
<b>Single Channel Models</b>				
P2	D10 DL10	DF10ST DFL10ST	<b>FP10001S</b>	0.2 - 2
P10	D10 DL10	DF10ST DFL10ST	<b>FP10002S</b>	1 - 10
P20	D200	DF30ST	<b>FP10003S</b>	2 - 20
P100	D200	DF100ST	<b>FP10004S</b>	10 - 100
P200	D200 D300	DF200ST DF300ST	<b>FP10005S</b>	20 - 200
P1000	D1000 D1200	DF1000ST DF1200ST	<b>FP10006S</b>	100 - 1000
P5000	D5000		<b>FP10007S</b>	500 - 5000
P10mL	D10mL		<b>FP10008S</b>	1000 - 10000
<b>Multichannel Models</b>				
P8x10	D10 DL10	DF10ST DFL10ST	<b>FP10009S</b>	0.5 - 10
P12x10			<b>FP10010S</b>	
P8x20	DL10 D200	DFL10ST DF30ST	<b>FP10011S</b>	2 - 20
P12x20			<b>FP10012S</b>	
P8x200	D200 D300	DF200ST DF300ST	<b>FP10013S</b>	20 - 200
P12x200			<b>FP10014S</b>	
P8x300	D200 D300	DF200ST DF300ST	<b>FP10015S</b>	20 - 300
P12x300			<b>FP10016S</b>	

Please contact your Gilson Account Manager for special text printing. Personalized pipettes will have the same part number with **SP** at the end.

[gilson.com/contactus](http://gilson.com/contactus)

LT800618/C | ©2022 Gilson S.A.S.  
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

