

MyPipetman®

SELECT

Precision That Adapts to Your Lab

MyPIPETMAN® Select combines precision control, ergonomic performance, and personalization, in one adaptable pipette that will be perfect for various protocols.



Set it. Lock it. Trust it.

With the Trilock, switch instantly between:

- **Free Setting Mode**
Quick, smooth adjustment when speed matters
- **Light Ratchet Mode**
Controlled precision
- **Locked Mode**
Secure volume locking for added confidence

Designed For All-Day Comfort

- Lightweight and well-balanced for all-day comfort
- Engineered with reduced pipetting forces
- Built to minimize hand strain and fatigue

Your Lab. Your Pipette.

- On-demand one-sided laser printing
- Dedicated pipettes to specific liquids or workflow steps to reduce cross-contamination risks
- Improved traceability and Lean lab organization



Ideal for Molecular biology, Microbiology, Cell Culture, Diagnostics, Pharmaceutical, Biochemistry, and General Laboratory Use

MyPIPETMAN® Select				
Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (µL)
Single Channel Models				
P2	D10 DL10	DF10ST DFL10ST	FP10001S	0.2 - 2
P10	D10 DL10	DF10ST DFL10ST	FP10002S	1 - 10
P20	D200	DF30ST	FP10003S	2 - 20
P100	D200	DF100ST	FP10004S	10 - 100
P200	D200 D300	DF200ST DF300ST	FP10005S	20 - 200
P1000	D1000 D1200	DF1000ST DF1200ST	FP10006S	100 - 1000
P5000	D5000		FP10007S	500 - 5000
P10mL	D10mL		FP10008S	1000 - 10000
Multichannel Models				
P8x10	D10 DL10	DF10ST DFL10ST	FP10009S	0.5 - 10
P12x10			FP10010S	
P8x20	DL10 D200	DFL10ST DF30ST	FP10011S	2 - 20
P12x20			FP10012S	
P8x200	D200 D300	DF200ST DF300ST	FP10013S	20 - 200
P12x200			FP10014S	
P8x300	D200 D300	DF200ST DF300ST	FP10015S	20 - 300
P12x300			FP10016S	

Please contact your Gilson Account Manager for special text printing. Personalized pipettes will have the same part number with SP at the end.

Ready-to-Use Pipette Kit

Start your experiment immediately with this complete solution including four pipettes, tips, and accessories. Everything you need in one box.

MyPIPETMAN® Select Kit				
Model	Part Number	Models	PIPETMAN® DIAMOND Tips	Accessories
MyPIPETMAN® Select 4-Pipette Kit	F167390	P2 P20 P200 P1000	DL10 D200 D1000	4 SINGLE pipette holders 4 Plastic ejectors 6 colored windows

Extra Accessories

Additional colored windows for customization and ease of use.

Colored Windows	
Model	Part Number
Window 3digit Dark Blue MyPIPETMAN	FP070101
Window 3digit White MyPIPETMAN	FP070102
Window 3digit Blue MyPIPETMAN	FP070103
Window 3digit Red MyPIPETMAN	FP070104
Window 3digit Yellow MyPIPETMAN	FP070105
Window 3digit Green MyPIPETMAN	FP070106

Discover MyPIPETMAN® Enterprise

Get Gilson quality with a solution tailored to your lab, MyPIPETMAN® Enterprise offers functional customization and personalization. Contact your Gilson Account Manager.



SCAN ME!

Learn more about MyPIPETMAN Select



Trilock™
3-position system for volume setting and locking



Free Volume Setting



Light Ratchet Volume Setting



Volume Locked

gilson.com/contactus

LT800618/E | ©2026 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

GILSON®