

# pipetman®

L

## Advantages Over Competition

### Ergonomics

- **Lightweight, well-balanced design:** Work longer with reduced fatigue and increased daily productivity at the bench
- **Low pipetting forces:** Minimizes strain and lowers the risk of repetitive strain injuries (RSIs)
- **Ambidextrous tip ejection:** Adjustable for left- or right-handed use, ensuring optimal comfort and flexibility

### Gilson Quality

- **Precision you can trust:** Delivers accurate and reproducible results, even in the most critical experiments
- **Built to last:** High-quality materials ensure long-term durability and reliable performance at the bench

### Safety

- **Secure volume locking system:** Prevents accidental adjustments and reduces the risk of costly errors
- **Fully autoclavable:** Protects against cross-contamination while preserving both pipette and user safety



#### Target Applications

- Biotechs
- Clinical labs
- High-intensity routine environments

# Battlecard

## PIPETMAN® L



PIPETMAN L Gilson
Single channel models + 0.2 µL - 10 mL
Multichannel models + 0.5 - 1200 µL
Fixed volume models + 1 µL to 5 mL

Volume-locking system
+ Lock-unlock

Tip Ejection System
+ Adjustable 3-position ejector, reduced forces

Autoclavable
+ Fully autoclavable without disassembly

Manufacturer Warranty
+ 3 years

RESEARCH PLUS Eppendorf	Tacta Sartorius	Pipet-Lite XLS+ Rainin MT	Finnpipette F1 Thermo	EVOLVE Integra	Transferpette S Brand	Research 3 Neo Eppendorf
Single channel models 0.1 µL - 10 mL	Single channel models 0.1 µL - 10 mL	Single channel models + 0.1 µL - 20 mL	Single channel models 0.2 µL - 10 mL	Single channel models 0.2 µL - 5 mL	Single channel models 0.1 µL - 10 mL	Single channel models 0.1 µL - 10 mL
Multichannel models + 0.5 µL - 1200 µL	Multichannel models 0.5 µL - 300 µL	Multichannel models + 0.5 µL- 1200 µL	Multichannel models 1 µL - 300 µL	Multichannel models 1 µL- 1250 µL	Multichannel models 0.5 µL - 300 µL	Multichannel models X No option
Fixed volume models 10 µL to 1000 µL	Fixed volume models X No option	Fixed volume models X No option	Fixed volume models + 10 µL to 10 mL	Fixed volume models X No option	Fixed volume models 10 µL to 2mL	Fixed volume models X No option

Volume-locking system						
X No system	+ Dual-function volume locking	+ Lever system with a brake (not absolute)	System with the push button	3-wheel adjustment with plunger down	Volume change protection with a slide button	+ Lock-unlock

Tip Ejection System						
Push button on the hook	Flat push button on the other side of the hook	Push button integrated in the handle	Push button on the other side of the hook	Push button above the counter	Push button on the other side of the hook	Push button on the other side of the hook

Autoclavable						
+ Fully autoclavable without disassembly	+ Fully autoclavable without disassembly	X Lower parts only (single channel only)	X Lower parts only (single channel only)	X Not recommended, lower parts only	+ Fully autoclavable without disassembly	+ Fully autoclavable without disassembly

Manufacturer Warranty						
3 years (local warranty)	2 years	+ 1 year	3 years (register pipettes for more)	+ 1 year	3 years single channel, 1 year multichannel	3 years